


Precision and Intuitive Photomanipulation by Mouse Click

- Intuitive system settings and precise experiment control using cellSens software
- Microsecond-accurate experiment synchronization and reproducibility
- Chromatically corrected diffraction-limited spot quality
- Full flexibility in position and shape of the bleaching or stimulation area
- Interactive “Click & Bleach” mode
- Integrated IX3 deck solution-easy to add to existing imaging platforms



cellFRAP – OLYMPUS Platform for Photomanipulation

Highly Accurate Control, Diffraction-Limited Stimulation, and the Ability to Precisely Reproduce Experiment Conditions

Understanding the complex dynamics inside cells by precise measurement of protein motions and interactions tells much about the local ongoing biological processes via, for example, the calculation of binding constants, coefficients of diffusion, and mobile fraction. Such analyses require a high level of accuracy and reproducibility. The Olympus cellFRAP system achieves such requirements and allows complex cellular processes to be investigated by photomanipulation. The system assures utmost local precision of manipulation through a diffraction-limited optical design. The microsecond-precise device control serves as a basis for best evaluation of fast processes and reproducibility of your experiments. Nicely adapted to the open source concept of the Olympus IX3 microscope, the cellFRAP can be combined with a variety of other advanced technologies—for example, TIRF or spinning-disc microscopy—and can easily be retrofitted to existing imaging platforms.



Standard Olympus cellFRAP-System setup with U-RTC



Advanced Olympus cellFRAP-System setup with Olympus 4L-TIRF system and U-RTCE

Precision and Reproducibility

Highest spatial and temporal resolution


The cellFRAP system guarantees highest spatial and temporal resolution for acquiring reliable and reproducible data in photomanipulation applications. Based on the system's diffraction-limited spot properties, bleaching areas can be controlled extremely accurately. Driven by the Olympus Real-time controller, the cellFRAP system offers an unrivaled short delay time of only 200 μ s between bleach and postbleach acquisition. Therefore, the system enables observation and collection of all important and immediate sample responses.

Diffraction-limited laser spot


The Olympus cellFRAP solution enables the best photomanipulation results without the need to move the sample. By utilizing a diffraction-limited laser spot, cellFRAP does not only guarantee the highest bleaching intensity but also ensures a precise and flexible choice in shape and position of the bleaching regions across the full field of view. High-quality chromatically corrected optics open up the possibility of using up to four different lasers during one experiment—for example, for optogenetics applications.

Reproducibility

Ultraprecise experimental control and reproducibility are of interest to most users. The cellFRAP system accounts for these demands by ultraprecise experimental control given by high system-level synchronization at ~ 1 μ s temporal resolution. In addition, the complete workflow of an experiment is documented and can easily be reproduced at any time.



Highest spatial resolution: A diffraction-limited spot size guarantees scanning precision and uniformity across the overall scan area



Highest temporal precision: 200 μ s between bleaching and acquisition

Ease of Use

Straight forward setup and control

Working with cellFRAP is intuitive and simple as the system is completely integrated into cellSens, the Olympus Life Science imaging platform. An automated calibration procedure makes it easy to perfectly align the scanner with the camera while taking into account the objectives and laser lines in use. Operating cellFRAP in the cellSens Graphical Experiment Manager offers intuitive system settings and experiment control, and it allows utmost operational flexibility and optimizes for the best speed performance. Thus, even complex acquisition workflows can be easily defined and reproduced, and any device operation can be automated.

Interactive experiment control

cellFRAP makes it possible to precisely adapt the bleaching regions to sample structures. An interactive definition of experimental conditions includes a free choice of bleaching/photomanipulation regions—variable in shape, including points and lines. Two different operation modes are available: Manipulation either based on predefined regions of interest or via interactive “Click & Bleach” mode. Signal intensity changes over time can also be observed online by displaying an intensity profile on demand.

Choose your method of evaluating data

In order to accurately evaluate your data, time markers for the start and end of the photomanipulation event are added to the recorded data. cellFRAP includes a qualitative FRAP analysis to classify your results. Here, recovery rates and mobile and immobile fractions can quickly be calculated from the recorded data including background and bleaching corrections. In addition, a Kymograph resembles a timeline visualization along a user-defined trajectory. An integrated export function allows for further evaluation of the data.


Expandable and Modular

Flexibility according to your needs

Integrated within the Olympus IX3 open-source concept, the cellFRAP deck solution can easily be adapted to an existing imaging platform and combines with a wide range of accessories like ZDC, climate chambers, manual or motorized stages, manipulators, etc. This allows for the setup of photomanipulation experiments in combination with other high-end applications—for example, spinning-disk and TIRF microscopy. The system can be equipped with up to four lasers for photomanipulation, a prerequisite for further sophisticated applications such as distinct photoswitching or optogenetics.




Just intuitive: Using the GEM (Graphical Experiment Manager) of cellSens, even complex acquisition and analysis workflows are defined and easily repeated




You can bleach points, lines or areas in variable shapes and modify your bleach conditions independently and separately for each shape



Qualitative FRAP analysis options are integrated into cellFRAP




Kymograph: Timeline visualization for easy interpretation of data and representation in publications





Example of a 'Click & Bleach' experiment

cellFRAP Dimensions and Applications



Lower port (deck 2)



Upper port (deck 1)

Applications

Photoactivation and photoswitching

Many cellular particles move quickly and unpredictably—for example, vesicles or intracellular pathogens such as viruses. For this reason, fluorescent molecules exhibiting photoactivation and photoconversion properties are used as markers to track the dynamic behavior. For example, fluorophores are available that can be changed-reversible or irreversible-in their spectral characteristics upon illumination with specific wavelengths—for example PA-GFP, KFP1, and Kaede. In particular, following the fate of, for example, a protein by an activated tagged molecule without the need to constantly providing additional fluorophores provides clear signals and minimal background conditions.

'Click & Bleach'

The 'Click & Bleach' function in the cellSens software enables the user to react to sample dynamics and to photomanipulate a target using direct mouse interaction. Multiple points can be bleached or activated and the resultant images can be captured during, or immediately after, the bleaching procedure.

FRAP

Fluorescence recovery after photobleaching (FRAP) provides an ideal method for calculating the coefficient of diffusion of a particular molecule within a certain target area of the cell. This is achieved


by photobleaching the fluorophore attached to a specific molecule within the target area and then assessing the dynamics of the return of fluorescence (recovery) to that area.

FLIP

Fluorescence loss in photobleaching (FLIP) is useful for studying molecular movement along cell membranes (lateral membrane fluidity) and membrane continuity, especially for membranous organelles. A defined area of the membrane or cell is continuously bleached, and the dynamic loss of fluorescence from the remaining area is measured based on the replacement of the bleached molecules in the defined area by unbleached ones from the remainder.

FLAP

Fluorescence localization after photobleaching (FLAP) is a clever adaptation of FRAP where the dynamics of a particular molecular species are of interest. The target molecule is coexpressed with two different fluorophores such that distribution and localization overlap. One of the fluorophores is then bleached at a defined location, and the movement of the molecule can be followed by comparing the relative and absolute distribution of the two fluorophores.



Schematic drawing of photo bleaching and advanced photomanipulation experiments



- OLYMPUS CORPORATION is ISO14001 certified.
- OLYMPUS CORPORATION is FM553994/ISO9001 certified.
- Illumination devices for microscope have suggested lifetimes. Periodic inspections are required. Please visit our website for details.
- This product is designed for use in industrial environments for the EMC performance. Using it in a residential environment may affect other equipment in the environment.
- All company and product names are registered trademarks and/or trademarks of their respective owners.
- Images on the PC monitors are simulated.
- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

www.olympus-lifescience.com

OLYMPUS

For enquiries - contact
www.olympus-lifescience.com/contact-us

OLYMPUS CORPORATION
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan
OLYMPUS EUROPA SE & CO. KG
Wendenstrasse 14-18, 20097 Hamburg, Germany
OLYMPUS SCIENTIFIC SOLUTIONS AMERICAS CORP.
48 Woerd Avenue, Waltham, MA 02453, U.S.A.
OLYMPUS SINGAPORE PTE LTD.
491B River Valley Road, #12-01/04 Valley Point Office Tower, Singapore 248373
OLYMPUS AUSTRALIA PTY. LTD.
3 Acacia Place, Notting Hill VIC 3168, Australia

OLYMPUS LATIN AMERICA, INC.
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126, U.S.A.
OLYMPUS (CHINA) CO., LTD.
A8F, Ping An International Financial Center, No. 1-3, Xinyuan South Road, Chaoyang District, Beijing, 100027 P.R.C.
OLYMPUS KOREA CO., LTD.
8F Olympus Tower, 446 Bongeunsa-ro, Gangnam-gu, Seoul, 135-509 Korea