

### OEM Microscope Components for Integration



# Contents

<b>1&gt;Welcome to UIS2 Optics</b> .....	1-1
<b>2.System Diagram</b>	
BX53M System Diagram (for Reflected and Reflected/Transmitted Light Combination) .....	2-1
BXFM System Diagram .....	2-3
BX63 System Diagram .....	2-5
BX53 System Diagram .....	2-7
BX43 System Diagram .....	2-9
BX3 Series Basic Motorized System Diagram .....	2-11
BX2 Series BXFM System Diagram .....	2-12
BX2 Series Motorized Unit System Diagram .....	2-13
BXFM-A System Diagram .....	2-13
<b>3.UIS2 Objectives</b>	
UIS2 Objectives (for Industrial Microscope) .....	3-1
M Plan Apochromat                   MPLAPON series .....	3-2
M Plan Apochromat                   MPLAPON100XO .....	3-2
M Plan SemiApochromat            MPLFLN series .....	3-3
Long WD M Plan SemiApochromat   LMPLFLN series .....	3-4
M Plan Achromat                    MPLN series .....	3-5
LCD Long WD M Plan SemiApochromat	
LCPLFLN-LCD series .....	3-6
Super Long WD M Plan Achromat    SLMPLN series .....	3-7
IR Long WD M Plan Achromat       LMPLN-IR series .....	3-8
IR M Plan Achromat                LCPLN-IR series .....	3-9
M Plan SemiApochromat BD        MPLFLN-BD series .....	3-10
M Plan SemiApochromat BDP       MPLFLN-BDP series .....	3-11
Long WD M Plan SemiApochromat BD	
LMPLFLN-BD series .....	3-12
M Plan Achromat BD               MPLN-BD series .....	3-13
White Light Interferometry Objective	
WL1100XMRTC .....	3-14
UIS2 Objectives (for Life Science Microscope) .....	3-15
Universal Plan Super Apochromat   UPLSAPO series .....	3-16
Plan Apochromat                   PLAPON series .....	3-18
Universal Plan Semi Apochromat/Plan Semi Apochromat	
UPLFLN, PLFLN series .....	3-19
Plan Achromat                    PLN series .....	3-21
Universal Plan Semi Apochromat for Phase Contrast	
UPLFLN-PH series .....	3-23
Plan Achromat for Phase Contrast   PLN-PH series .....	3-24
Universal Plan Semi Apochromat for Polarizing	
UPLFLN-P series .....	3-25
Achromat for Polarizing            PLN-P, ACHN-P series .....	3-26
Plan Achromat (ND)                PLN-CY, PLFLN-CY series ...	3-27
Long Working Distance Universal Plan Semi Apochromat	
LUCPLFLN series .....	3-28
Long Working Distance Universal Plan Semi Apochromat for Relief Contrast	
CPLFLN-RC, LUCPLFLN-RC series .....	3-29
Long Working Distance Universal Plan Semi Apochromat for Phase Contrast	
CPLFLN-PH, LUCPLFLN-PH series .....	3-30
Culture Specimen Objectives for Phase Contrast	
CPLN-PH, LCACHN-PH series .....	3-31
Culture Specimen Objectives for Relief Contrast	
CPLN-RC, LCACHN-RC series .....	3-32
No Cover Water Immersion for Fixed Stage Upright Microscope	
UMPLFLN-W, LUMPLFLN-W series .....	3-33
No Cover Water Immersion for Fixed Stage Upright Microscope	
XLUMPLFLN20XW .....	3-34
Universal Apochromat               UAPON 340 series .....	3-35
TIRF Objectives                    APON, UAPON series .....	3-36

<b>4.Microscope Frame</b>	
BX53M: Upright Transmitted & Reflected Light Microscope Frame	
BX53MTRF-S .....	4-1
BX53M: Upright Reflected Light Microscope Frame	
BX53MRF-S .....	4-2
BX3: Automated Transmitted Light Microscope Frame	
BX63F .....	4-3
BX3: Semi-Motorized Fluorescence Transmitted Light Microscope Frame	
BX53F .....	4-4
BX3: Manual System Transmitted Light Microscope Frame	
BX43F .....	4-5
BX3: Transmitted Ergonomic Microscope Frame	
BX46F .....	4-6
BX2: Upright Motorized Transmitted/Reflected Frame	
BX61TRF .....	4-7
BXFM Frame                        BXFM-F .....	4-8
BXFM System Configuration Example1	
BXFM-F + BXFM-ILH + BXFM-ILHSPU ....	4-9
BXFM System Configuration Example2	
BXFM-F + BXFM-ILHS .....	4-10
Stands for BXFM .....	4-11
<b>5.Illumination Units</b>	
Reflected Light Illuminator for BX53M .....	5-1
Coded Reflected Light Illuminator for BX53M Frame .....	5-2
Reflected Illuminator for BX3 Series .....	5-3
Reflected Light Illuminator for BX2 series .....	5-4
Mounting Dimensions of Illuminator (BX3M-RLA-S, BX3M-RLAS-S, BX3M-URAS-S, BX3M-KMA-S, BX3-RFAS, BX3-URA)	
Mounting Dimensions of Illuminators (BX-RLA2 and BX-URA2) .....	5-5
Compact Reflected Light Illuminator for BF .....	5-6
<b>6.Light Source Units</b>	
LED Lamp Housing for BX53M .....	6-1
LED Lamp Housing for BX3 Series .....	6-2
Lamp Housings .....	6-3
Halogen Illumination .....	6-4
Halogen Fiber Illumination Accessories .....	6-5
Lamp Housing Accessory .....	6-6
<b>7.Condenser Units</b>	
Universal Condenser .....	7-1
Condenser .....	7-2
<b>8.Observation Tubes</b>	
Super Widefield Trinocular Observation Tubes .....	8-1
Widefield Trinocular Observation Tubes .....	8-2
Single Port Tube with Lens .....	8-2
Tilting Binocular Tube .....	8-3
Binocular Tube .....	8-4
<b>9.Intermediate Tubes Accessories</b>	
Intermediate Tubes .....	9-1
Dual port tube with C mounts .....	9-3
<b>10.Eyepieces</b>	
Eyepieces .....	10-1
<b>11.Revolving Nosepieces</b>	
Revolving Nosepieces for BF Objectives .....	11-1
Revolving Nosepieces for BF/DF Objectives .....	11-2
Coded Sextuple Revolving Nosepiece .....	11-3
<b>12.Video Camera Adapters</b>	
C-mount Video Camera Ports .....	12-1
Video Camera Mount Adapters .....	12-2
Video Camera Port .....	12-2
<b>13.Motorized Units</b>	
Motorized Universal Reflected Illuminator for BX2 Series .....	13-1
Motorized Units .....	13-2
Control Box for BX2 Series .....	13-3
Motorized Units for BX2 Series .....	13-4
Motorized Modular Microscope .....	13-5
Motorized Units for BX3 Series .....	13-6
Control Box for BX3 Series .....	13-7
Control Box for BX53M/BXFM .....	13-8
<b>14.Optical Terminology</b> .....	14-1


# Welcome to UIS2 Optics

## UIS2:

### The System That Maximizes the Advantage of Infinity-corrected Optics

#### What is Infinity-corrected Optics?

UIS2 optics is an infinity-corrected optical system — in other words, a system in which light passes from the specimen through the objectives without forming an image along the way. Instead, it travels in the form of infinity parallel rays to the tube lens. The tube lens is where the intermediate image is formed, whereas in finite-corrected optics, this is done by the objective.



#### Advantages of Infinity-corrected Optics


This system, known as "infinity-corrected optics", offers a number of advantages:

- There is no change in magnification even when the distance between the objective and tube lens is altered.
- With the total magnification remaining constant, there is no image aberration — even when prisms or sliders are interposed between the objectives and the tube lens.

As thousands of users have found by experience, these advantages are crucial to composing the ideal microscope optical system. What's more, it is even possible to freely insert or remove intermediate attachments in the parallel rays of light between the objectives and tube lens, allowing the creation of

user-specific or task-specific optical systems. To establish real flexibility with such a system, it is necessary to eliminate the occurrence of coma aberration.


*\*In UIS2 Objectives, the parfocal distance is designed at 45 mm and the focal length of the tube lens is 180 mm.*



#### Basic Dimensions of UIS2 Optical System

The UIS2 optical system optimally corrects aberration with a dedicated telan lens and an eyepiece so that the coma aberration and flatness are not degraded even when the telan lens exit pupil position is changed by changing the objective and telan distance. This makes it possible to use a distance of 50 mm to 170 mm from objective mounting position to the single port tube with lens.

*\*Coma aberration: refer to the optical terminology at the end of this document.*



## Features of UIS2 Objectives

UIS2 objective lenses provide compatibility (screw diameter, optical performance) with the UIS optical system and have the following features compared to conventional objectives.

### **1. Wavefront Aberration Control**

The Olympus UIS2 objectives set a new standard, with wavefront aberration control in addition to common performance standards of NA and W.D. Olympus challenges farther advanced optics which has not been fulfilled by the conventional standards. We offer splendid performance objectives by minimizing the aberrations that lower resolution.

*\*Wave front aberration: refer to the optical terminology at the end of this document.*

### **2. Objective Lenses with Splendid Image Parcentricity**

High power SemiAchromatic UIS2 objectives make the centration tolerance between objectives on the microscope nosepiece keep the image within the center of the field of view even with digital cameras. (50X or higher power in both MPLFLN and LMPLFLN series)

### **3. Improvement of Color Reproducibility**

UIS2 objectives realize true color reproduction without any chromatic shifts using stringently selected high transmittance glass and advanced coating technology that provides high transmittance which is flat over a wide band wavelength. In addition, since the total optical system, including the tube lens is designed to reproduce a true color, clear images faithful to the specimen are obtained even with digital imaging.

### **4. Lightening**

Weight has been reduced to approximately 2/3 that of conventional products by using an aluminum objectives barrel cover. This has the effect of lightening the load on the devices at objective up/down, suppressing vibrations by lowering the inertial force at objective switching, etc. (MPLFLN series, LMPLFLN series)

### **5. Adoption of Eco-lens**


The glass materials of UIS2 objectives are all lead- and cadmium-free eco-glass.

Based on our conviction that the UIS2 system maximizes the advantages of infinity-corrected optical systems, we confidently recommend the UIS2-featured Olympus microscope units for all your high-precision needs in research, inspection and production equipment.

*\*Refer to the Olympus home page for detailed objective lenses specifications.*


# System Diagram

## BX53M System Diagram (for Reflected and Reflected/Transmitted Light Combination)




\*For transmitted light combination only


## Control box and cable connection diagram



## Stand-alone light manager configuration




## MIX observation configuration




# System Diagram

## BXFM System Diagram




\*U-ST is not available with BX3M-ILH.


### Control box and cable connection diagram




### MIX observation configuration



### Motorized nosepiece configuration









# System Diagram

## BX53 System Diagram




\*1 Slight vignetting may occur in combination with an additional intermediate attachment or observation method. \*2 Require an additional intermediate attachment or fluorescence illuminator. \*3 Cannot be used with U-TTLBI. \*4 Compatible with FN 22. \*5 Cannot be used with BX3-URA. \*6 Stand is a standard equipment of the U-MDOSV and U-MDO10R3. \*7 An auxiliary lens is equipped.




# System Diagram

## BX43 System Diagram






\*1 Slight vignetting may occur in combination with an additional intermediate attachment or observation method. \*2 Require an additional intermediate attachment or fluorescence illuminator.

\*3 Cannot be used with U-TTLBI. \*4 Compatible with FN 22. \*5 An auxiliary lens is equipped.



# System Diagram


## BX3 Series Basic Motorized System Diagram



\*Please refer to "Section xx Motorized Unit" for each motorized unit in the detail.

\*\*Please consult your nearest Olympus representative in your region about the detailed motorized system configuration and combination.


## BX2 Series BAXM System Diagram




\*Different types may be offered in each area.

# System Diagram


## BX2 Series Motorized Unit System Diagram




## BXFM-A System Diagram



## Definition for Objective Abbreviations



## Objective Notation



## Objective series List

Series	Magnification	BF	DF	DIC*1	POL	FL	FN	Remarks
MPLAPON	50/100	○		○ U	○		26.5	
MPLAPON O	100	○					26.5	
MPLFLN	1.25/2.5	○					1.25X: 22/2.5X: 26.5	Use together with polarizer and analyzer recommended.
	5/10/20/40*2/50/100	○		○ U	○	○*3	26.5	
LMPLFLN	5/10/20/50/100	○		○ L	○	○	26.5	
MPLN	5/10/20/50/100	○					22	
LCPLFLN-LCD	20/50/100	○		○ L			26.5	for LCD
SLMPLN	20/50/100	○					26.5	
LMPLN-IR	5/10	○					22	for near-IR observation
LCPLN-IR	20/50/100	○					22	for near-IR observation
MPLFLN-BD	5/10/20/50/100/150	○	○	○ U	○	○*3	26.5	
MPLFLN-BDP	5/10/20/50/100	○	○	○ U	○	○*3	26.5	
LMPLFLN-BD	5/10/20/50/100	○	○	○ L	○	○	26.5	
MPLN-BD	5/10/20/50/100	○	○				22	
WL100XMRTC	100X						22	Mirau objective

\*1 DIC Prism U-DICR: UM/LM Position, U-DICRHC: LM Position Fixed, U-DICRH: UM Position Fixed. \*2 40X: BF Only \*3 5-20X: U Excitation Also Possible  
 ○: Responds ◎: Optimally Responds BF: Brightfield DF: Darkfield DIC: Differential Interference Contrast POL: Polarized light FL: Fluorescence

## Features of Objective Series

### ● MPLAPON series: M Plan Apochromat — P 3-2

Plan Apochromat objectives that correct chromatic aberrations at optimal level among Olympus objectives. Olympus guarantees\* optical performance (wavefront aberration) with a Strehl ratio\*\* of 95% or better. They are also designed for use with Olympus' U-AFA2M active AF unit.

### ● MPLAPON100XO: M Plan Apochromat — P 3-2

This is a Plan Apochromat objective of the oil-immersion type\*\*\* that features a numerical aperture of 1.4. It provides excellent chromatic aberration correction and a high resolving power.

### ● MPLFLN series: M Plan SemiApochromat — P 3-3

Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. The lineup consists of 8 objectives ranging from 1.25X to 100X, and provides a W.D. of 1 mm or longer (except 40X). Since the exit pupil position of the 5X-100X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification (40X is not applicable to DIC observation). For very low magnifications (1.25X, 2.5X), use together with analyzer and polarizer of the reflected light illuminator.

### ● LMPLFLN series: Long WD M Plan SemiApochromat — P 3-4

Long working distance Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. Suitable with samples having a height difference and in preventing collision, as the working distance is long. Also, since the exit pupil position of the 5X-100X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification.

### ● MPLN series: M Plan Achromat — P 3-5

Plan Achromat objectives providing excellent image flatness up to FN 22.

### ● LCPLFLN-LCD series: LCD Long WD M Plan SemiApochromat — P 3-6

Differentiated objective series for observation of LCD panels and other samples through a glass substrate. Aberration correction matched to the glass thickness is accomplished using a correction collar.

### ● SLMPLN series: Super Long WD M Plan Achromat — P 3-7

Plan Achromat objectives with high magnification and super long working distance. Three magnifications, 20X, 50X and 100X are available. For 5X or 10X objectives, select from the LMPLFLN series.

### ● LMPLN-IR series: IR Long WD M Plan Achromat — P 3-8

Objective series exclusive for the near-infrared microscopy largely of the internal structure in silicon wafers.

### ● LCPLN-IR series: IR Long WD M Plan Achromat — P 3-9

Exclusive for the near-infrared microscopy largely of the internal structure in silicon wafers, LCPLN-IR series has correction collar for aberration dependent on thickness of silicon or glass substrate.

### ● MPLFLN-BD series: M Plan SemiApochromat BD — P 3-10

Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. The series provides a W.D. of 1 mm or longer. Since the exit pupil position of the 5X-150X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification.

### ● MPLFLN-BDP series: M Plan SemiApochromat BDP — P 3-11

Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. The series provides a W.D. of 1 mm or longer. Since the exit pupil position of the 5X-100X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification. The BDP series optimizing brightfield/darkfield and polarized light characteristics is appropriate for Nomarski DIC and polarized light observations.

### ● LMPLFLN-BD series: Long WD M Plan SemiApochromat BD — P 3-12

Long working distance Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. Suitable with samples having a height difference and in preventing collision, as the working distance is long. Also, since the exit pupil position of the 5X-100X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification.

### ● MPLN-BD series: M Plan Achromat BD — P 3-13

Plan Achromat objectives providing excellent image flatness up to FN 22.

### ● WL100XMRTC series: White Light Interferometry Objective — P 3-14

This objective is designed for the Mirau style of white light interferometers and maintains a high level of temperature tolerance. The optimized NA of 0.8 provides improved light gathering, with a working distance of 0.7 mm.

*\*Definition of guaranteed values: Measurements assessed with Olympus' Interferometer for Transmitted Wavefront Measurement under specified conditions (measurement: temp. 23°C±1°C; assessment: measurement within the 97% range of the pupil dia.).*

*\*\*Strehl ratio: Indicates in percent (%) the ratio of the proportion of light that an actual optical system can concentrate with respect to the proportion of light concentrated in the image plane (central intensity) by an ideal, aberration-free optical system, with the latter serving as 100%. A higher percentage indicates a higher quality optical system.*


*\*\*\*Specified oil: IMMOIL-F30CC*

**M Plan Apochromat**


**MPLAPON series**

Plan Apochromat objectives that correct chromatic aberrations at optimal level among Olympus objectives. Olympus guarantees\* optical performance (wavefront aberration) with a Strehl ratio\*\* of 95% or better. They are also designed for use with Olympus' U-AFA2M active AF unit.

**MPLAPON50X**



**MPLAPON100X**



\* Definition of guaranteed values: Measurements assessed with Olympus' Interferometer for Transmitted Wavefront Measurement under specified conditions (measurement: temp. 23°C±1°C; assessment: measurement within the 97% range of the pupil dia.).

\*\* Strehl ratio: Indicates in percent (%) the ratio of the proportion of light that an actual optical system can concentrate with respect to the proportion of light concentrated in the image plane (central intensity) by an ideal, aberration-free optical system, with the latter serving as 100%. A higher percentage indicates a higher quality optical system.

Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super Widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)
MPLAPON50X	0.95	0.35	3.6	139	500	0.44	1	500	0.53	1
MPLAPON100X	0.95	0.35	1.8	125	1000	0.22	0.67	1000	0.27	0.7


Screw: W20.32 × 0.706 (0.8" × 1/36")

**M Plan Apochromat**

**MPLAPON100XO**

This is a Plan Apochromat objective of the oil-immersion type\*\*\* that features a numerical aperture of 1.4. It provides excellent chromatic aberration correction and a high resolving power.

**MPLAPON100XO**



\*\*\* Specified Oil: IMMOIL-F30CC

Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super Widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)
MPLAPON100XO	1.4	0.1	1.8	158	1000	0.22	0.59	1000	0.27	0.59


Screw: W20.32 × 0.706 (0.8" × 1/36")

**M Plan SemiApochromat**


**MPLFLN series**

Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. The lineup consists of 8 objectives ranging from 1.25X to 100X, and provides a W.D. of 1 mm or longer (except 40X). Since the exit pupil position of the 5X–100X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification (40X is not applicable to DIC observation). For very low magnifications (1.25X, 2.5X), use together with analyzer and polarizer of the reflected light illuminator.


**MPLFLN1.25X**




**MPLFLN2.5X**




**MPLFLN5X**




**MPLFLN10X**




**MPLFLN20X**




**MPLFLN40X**



**MPLFLN50X**



**MPLFLN100X**




Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super Widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)
MPLFLN1.25X	0.04	3.5	145	122	12.5	17.6	870	—	—	—
MPLFLN2.5X	0.08	10.7	72	106	25	8.8	220	25	10.6	220
MPLFLN5X	0.15	20	36	51.5	50	4.4	59	50	5.3	59
MPLFLN10X	0.3	11	18	68.1	100	2.2	15	100	2.7	15
MPLFLN20X	0.45	3.1	9	70.4	200	1.1	5.2	200	1.3	5.1
MPLFLN40X	0.75	0.63	4.5	120	400	0.55	1.66	400	0.66	1.66
MPLFLN50X	0.8	1	3.6	89.9	500	0.44	1.3	500	0.53	1.3
MPLFLN100X	0.9	1	1.8	90.9	1000	0.22	0.73	1000	0.27	0.73


Screw: W20.32 × 0.706 (0.8" × 1/36")

Long working distance Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. Suitable with samples having a height difference and in preventing collision, as the working distance is long. Also, since the exit pupil position of the 5X–100X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification.


**LMPLFLN5X**




**LMPLFLN10X**




**LMPLFLN20X**



**LMPLFLN50X**



**LMPLFLN100X**




Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super Widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (µm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (µm)
LMPLFLN5X	0.13	22.5	36	50	50	4.4	70	50	5.3	70
LMPLFLN10X	0.25	21	18	54	100	2.2	18	100	2.7	18
LMPLFLN20X	0.4	12	9	73	200	1.1	6.1	200	1.3	6.1
LMPLFLN50X	0.5	10.6	3.6	77	500	0.44	2.5	500	0.53	2.5
LMPLFLN100X	0.8	3.4	1.8	94	1000	0.22	0.87	1000	0.27	0.87


Screw: W20.32 × 0.706 (0.8" × 1/36")

Plan Achromat objectives providing excellent image flatness up to FN 22.


**MPLN5X**




**MPLN10X**




**MPLN20X**



**MPLN50X**



**MPLN100X**



Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus ( $\mu$ m)
MPLN5X	0.1	20	36	64	50	4.4	98
MPLN10X	0.25	10.6	18	80	100	2.2	18
MPLN20X	0.4	1.3	9	111	200	1.1	6.1
MPLN50X	0.75	0.38	3.6	13	500	0.44	1.4
MPLN100X	0.9	0.21	1.8	116	1000	0.22	0.73


Screw: W20.32 x 0.706 (0.8" x 1/36")

# LCD Long WD M Plan SemiApochromat


## LCPLFLN-LCD series

Differentiated objective series for observation of LCD panels and other samples through a glass substrate. Aberration correction matched to the glass thickness is accomplished using a correction collar.


**LCPLFLN20XLCD**



**LCPLFLN50XLCD**



**LCPLFLN100XLCD**



\*Value at glass thickness 0.7 mm observation.

Unit: mm

Objective	LCPLFLN20XLCD			LCPLFLN50XLCD			LCPLFLN100XLCD		
Corresponding Glass Thickness (mm)	0-1.2			0-1.2			0-0.7		
Correction Collar Indication	0	0.7	1.2	0	0.7	1.2	0	0.5	0.7
W.D. (mm)	8.3	7.8	7.4	3	2.5	2.2	1.2	0.98	0.9
Correction System	Correction Collar			Correction Collar			Correction Collar		

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D.** (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)
LCPLFLN20XLCD	0.45	7.8	9	146	200	1.1	5.2	200	1.3	5.2
LCPLFLN50XLCD	0.7	2.5	3.6	170	500	0.44	1.6	500	0.53	1.6
LCPLFLN100XLCD	0.85	0.9	1.8	185	1000	0.22	0.79	1000	0.27	0.79


Screw: W20.32 × 0.706 (0.8" × 1/36")

\*\*The figure shown here is the value when the correction collar indication is 0.7.

## Super Long WD M Plan Achromat

### SLMPLN series

Plan Achromat objectives with high magnification and super long working distance. Three magnifications, 20X, 50X and 100X are available. For 5X or 10X objectives, select from the LMPLFLN series.




Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super Widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)
SLMPLN20X	0.25	25	9	56	200	1.1	11.4	200	1.3	11.4
SLMPLN50X	0.35	18	3.6	74	500	0.44	4.2	500	0.53	4.2
SLMPLN100X	0.6	7.6	1.8	100	1000	0.22	1.3	1000	0.27	1.3


Screw: W20.32 × 0.706 (0.8" × 1/36")

Objective series exclusive for the near-infrared microscopy largely of the internal structure in silicon wafers.

**LMPLN5XIR**



**LMPLN10XIR**




Unit: mm

UIS Objectives				
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)
LMPLN5XIR	0.1	23	36	55
LMPLN10XIR	0.3	18	18	78

Screw: W20.32 × 0.706 (0.8" × 1/36")

Exclusive for the near-infrared microscopy largely of the internal structure in silicon wafers, LCPLN-IR series has correction collar for aberration dependent on thickness of silicon or glass substrate.



Unit: mm

**Silicon thickness correction**

Objective	LCPLN20XIR			LCPLN50XIR			LCPLN100XIR		
Corresponding Silicon Thickness (mm)	0-1.2			0-1.2			0-1.0		
Correction Collar Indication	0	0.7	1.2	0	0.6	1.2	0	0.5	1
W.D.* (mm)	8.3	8.2	8	4.5	4.3	4.1	1.2	1.1	1
Correction System	Correction Collar			Correction Collar			Correction Collar		

\*With the use of 1100 nm laser.

**Silicon thickness correction**

Objective	LCPLN20XLCD			LCPLN50XLCD		LCPLN100XLCD	
Corresponding Glass Thickness (mm)	0-1.2			0-1.2		0-0.7	
Correction Collar Indication	0	0.7	1.2	0	1.2	0	0.7
W.D.* (mm)	8.3	7.9	7.6	4.5	3.7	1.2	0.9
Correction System	Correction Collar			Correction Collar		Correction Collar	

\*With the use of 1064 nm laser.


UIS2 Objectives				
Objective (magnification)	NA*	W.D.* (mm)	Focal Distance f (mm)	Weight (g)
LCPLN20XIR	0.45	8.3	9	149
LCPLN50XIR	0.65	4.5	3.6	169
LCPLN100XIR	0.85	1.2	1.8	184

Screw: W20.32 × 0.706 (0.8" × 1/36")


\*The figure shown here is the value when the correction collar indication is 0.

Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. The series provides a W.D. of 1 mm or longer. Since the exit pupil position of the 5X-150X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification.


**MPLFLN5XBD**




**MPLFLN10XBD**




**MPLFLN20XBD**




**MPLFLN50XBD**



**MPLFLN100XBD**



**MPLFLN150XBD**



Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super Widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)
MPLFLN5XBD	0.15	12	36	95.5	50	4.4	59	50	5.3	59
MPLFLN10XBD	0.3	6.5	18	82.8	100	2.2	15	100	2.7	15
MPLFLN20XBD	0.45	3	9	87.7	200	1.1	5.2	200	1.3	5.2
MPLFLN50XBD	0.8	1	3.6	99.8	500	0.44	1.3	500	0.53	1.3
MPLFLN100XBD	0.9	1	1.8	98.9	1000	0.22	0.73	1000	0.27	0.73
MPLFLN150XBD	0.9	1	1.2	104.8	1500	0.15	0.6	1500	0.18	0.6


Screw: W26 × 0.706

**M Plan SemiApochromat BDP**


**MPLFLN-BDP series**

Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. The series provides a W.D. of 1 mm or longer. Since the exit pupil position of the 5X-100X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification. The BDP series optimizing brightfield/darkfield and polarized light characteristics is appropriate for Nomarski DIC and polarized light observations.


**MPLFLN5XBDP**




**MPLFLN10XBDP**




**MPLFLN20XBDP**



**MPLFLN50XBDP**



**MPLFLN100XBDP**




Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super Widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)
MPLFLN5XBDP	0.15	12	36	95.5	50	4.4	59	50	5.3	59
MPLFLN10XBDP	0.25	6.5	18	83.3	100	2.2	18	100	2.7	18
MPLFLN20XBDP	0.4	3	9	88.5	200	1.1	6.1	200	1.3	6.1
MPLFLN50XBDP	0.75	1	3.6	100.5	500	0.44	1.4	500	0.53	1.4
MPLFLN100XBDP	0.9	1	1.8	101.5	1000	0.22	0.73	1000	0.27	0.73


Screw: W26 × 0.706

Long working distance Plan SemiApochromat objectives, giving high-level correction for chromatic aberration. Suitable with samples having a height difference and in preventing collision, as the working distance is long. Also, since the exit pupil position of the 5X-100X objectives is standardized, the position of the DIC prism does not have to be switched when changing the magnification.


**LMPLFLN5XBD**




**LMPLFLN10XBD**




**LMPLFLN20XBD**



**LMPLFLN50XBD**



**LMPLFLN100XBD**




Unit: mm

UIS2 Objectives					Widefield Eyepiece WHN10X FN 22			Super Widefield Eyepiece SWH10X FN 26.5		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (μm)
LMPLFLN5XBD	0.13	15	36	81	50	4.4	70	50	5.3	70
LMPLFLN10XBD	0.25	10	18	84	100	2.2	18	100	2.7	18
LMPLFLN20XBD	0.4	12	9	86	200	1.1	6.1	200	1.3	6.1
LMPLFLN50XBD	0.5	10.6	3.6	85	500	0.44	2.5	500	0.53	2.5
LMPLFLN100XBD	0.8	3.3	1.8	102	1000	0.22	0.87	1000	0.27	0.87


Screw: W26 × 0.706

Plan Achromat objectives providing excellent image flatness up to FN 22.


**MPLN5XBD**




**MPLN10XBD**




**MPLN20XBD**



**MPLN50XBD**



**MPLN100XBD**



Unit: mm


UIS2 Objectives					Widefield Eyepiece WHN10X FN 22		
Objective (magnification)	NA	W.D. (mm)	Focal Distance f (mm)	Weight (g)	Total Magnifications	Practical Field of View (mm)	Depth of Focus (µm)
MPLN5XBD	0.1	12	36	137	50	4.4	98
MPLN10XBD	0.25	6.5	18	155	100	2.2	18
MPLN20XBD	0.4	1.3	9	162	200	1.1	6.1
MPLN50XBD	0.75	0.38	3.6	157	500	0.44	1.4
MPLN100XBD	0.9	0.21	1.8	160	1000	0.22	0.73

Screw: W26 × 0.706

**White Light Interferometry Objective**

**WLI100XMRTC**

This objective is designed for the Mirau style of white light interferometers and maintains a high level of temperature tolerance. The optimized NA of 0.8 provides improved light gathering, with a working distance of 0.7 mm.




unit: mm


Objective (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
WLI100XMRTC	0.8	0.7	0	—	—	—	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

## Definition for Objective Abbreviations



## Objective Notation



## Objective series List

### Objective series for standard biological samples

Series	Magnification	BF	DF	DIC*	POL	FL	FN	Remarks
UPLSAPO	4X/10X/20X/20XO/40X/60XO/60XW/100XO	○	10X/20X/20XO	○ (except 4X)	○	○	26.5	
PLAPON	1.25X/2X/60XO	○		60XO	60XO	60XO	26.5	
UPLFLN	4X/10X/20X/40X/40XO/60X/60XOI/100XO/100XOI	○	10X/20X/40X/60XOI/100XOI	○ (except 4X)	○	○	26.5	
PLFLN	100X	○				○	26.5	
PLN	2X/4X/10X/20X/40X/50XOI/100XO	○	10X/20X/40X/50XOI	10X/20X/40X/50XOI		○	22	
UPLFLN-PH	4XPH/10XPH/20XPH/40XPH/60XOIPH/100XOPH	○	10XPH/20XPH/40XPH/60XOIPH			○	26.5	
PLN-PH	10XPH/20XPH/40XPH/100XOPH	○	○ (except 100XOPH)				22	
UPLFLN-P	4XP/10XP/20XP/40XP/100XOP	○	10XP/20XP/40XP	○ (except 4XP)	○	○	26.5	
PLN-P/ACHN-P	4XP/10XP/20XP/40XP/100XOP	○	10XP/20XP/40XP		○	○	22	
PLN-CY	2XCX/4XCX/10XCX/20XCX/(FLN) 10XCX	○					22	PLFLN10XCX FN 26.5

\*These objectives are suitable for standard biological samples embedded on slide glass with 0.17 mm cover slip and mainly used with upright microscopes.

### Objective series for cultured samples

Series	Magnification	BF	DF	DIC	POL	FL	FN	Remarks
LUCPLFLN	20X/40X/60X	○	○	○	○	○	22	
LUCPLFLN-RC/ UCPLFLN-RC	10XRC/20XRC/40XRC	○	○			○	22	
LUCPLFLN-PH/ UCPLFLN-PH	10XPH/20XPH/40XPH/60XPH	○	○			○	22	
CPLN-PH/ LCACHN-PH	10XPH/20XPH/40XPH	○	○				22	
CPLN-RC/ LCACHN-RC	10XRC/20XRC/40XRC	○	○				22	

These objectives are suitable for cultured tissue/cell observation in dish/ bottle or micro-plate and mainly used with inverted microscopes.

### Objective series for special applications


Series	Magnification	BF	DF	DIC	POL	FL	FN	Remarks
LUMPLFLN/ UMPLFLN	10XW/20XW/40XW/60XW	○	10XW/20XW	○	○	○	26.5	
XLUMPLFLN	20XW	○		○		○	22	Perofocal length 75 mm
APON 340	20XW/40XO/40XW	○	20XW/40XW	○	○	○	22	
TIRF Objectives	60XO/100XHO/100XO/150XO	○		○	○	○	22	

## Features of Objective Series (Please refer to the following page in the detail for each Objective.)


- UPLSAPO: Universal Plan Super Apochromat — P 3-16
- PLAPON: Plan Apochromat — P 3-18
- UPLFLN: Universal Plan Semi Apochromat/Plan Semi Apochromat — P 3-19
- PLN: Plan Achromat — P 3-21
- UPLFLN-PH UPlanFI-P Universal Plan Semi Apochromat for Phase Contrast — P 3-23
- PLN-PH: Plan Achromat for Phase Contrast — P 3-24
- UPLFLN-P: UPlanFI-P Universal Plan Semi Apochromat for Polarizing. — P 3-25
- PLN/ACHN-P: Achromat for Polarizing. — P 3-26
- PLN-CY: Plan Achromat (ND) — P 3-27
- LUCPLFLN, UCPLFLN: Long Working Distance Universal Plan Semi Apochromat — P 3-28
- LUCPLFLN/UCPLFLN-RC: Long Working Distance Universal Plan Semi Apochromat for Relief Contrast — P 3-29
- LUCPLFLN/UCPLFLN-PH: Long Working Distance Universal Plan Semi Apochromat for Phase Contrast — P 3-30
- CPLN/LCACHN-PH: Culture Specimen Objectives for Phase Contrast — P 3-31
- CPLN/LCACHN-RC: Culture Specimen Objectives for Relief Contrast — P 3-32
- LUMPLFLN/UMPLFLN: No Cover Water Immersion for Fixed Stage Upright Microscope — P 3-33
- XLUMPLFLN: No Cover Water Immersion for Fixed Stage Upright Microscope — P 3-34
- APON 340: Universal Apochromat — P 3-35
- TIRF Objectives — P 3-36

The UPLSAPO Super Apochromat objectives fully compensate for both spherical and chromatic aberrations from the UV to the near-infrared region. The high sensitivity to fluorescence emissions allows the acquisition of sharp, clear images, without color shift, even in brightfield and Nomarski DIC observations. For quality and performance, they offer an unbeatable solution to every kind of digital imaging need.


**UPLSAPO4X**




**UPLSAPO10X2**



**UPLSAPO20X**



**UPLSAPO20XO**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
UPLSAPO4X	0.16	13.0	—	—	—	U/BG/IR	26.5
UPLSAPO10X2	0.40	3.1	0.17	—	—	U/BG/IR	26.5
UPLSAPO20X	0.75	0.6	0.17	—	Yes	U/BG/IR	26.5
UPLSAPO20XO	0.85	0.2	—	Oil	Yes	U/BG/IR	26.5

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.

## Universal Plan Super Apochromat


### UPLSAPO series

The UPLSAPO Super Apochromat objectives fully compensate for both spherical and chromatic aberrations from the UV to the near-infrared region. The high sensitivity to fluorescence emissions allows the acquisition of sharp, clear images, without color shift, even in brightfield and Nomarski DIC observations. For quality and performance, they offer an unbeatable solution to every kind of digital imaging need.


**UPLSAPO40X2**




**UPLSAPO60XO**



**UPLSAPO60XW**



**UPLSAPO100XO**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
UPLSAPO40X2	0.95	0.18	0.11-0.23	—	Yes	U/BG/IR	26.5
UPLSAPO60XO	1.35	0.15	0.17	Oil	Yes	U/BG/IR	26.5
UPLSAPO60XW	1.20	0.28	0.13-0.21	Water	Yes	U/BG/IR	26.5
UPLSAPO100XO	1.40	0.13	0.17	Oil	Yes	U/BG/IR	26.5

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.

**Plan Apochromat**


**PLAPON series**

The PLAPON Apochromat objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. The high sensitivity to fluorescence emissions allows the acquisition of sharp, clear images, without color shift, even in brightfield and Nomarski DIC observations. For quality and performance, they offer an unbeatable solution to every kind of digital imaging need.


**PLAPON1.25X**



**PLAPON2X**



**PLAPON60XO**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
PLAPON1.25X	0.04	5.0	—	—	—	(BG)/(IR)	26.5
PLAPON2X	0.08	6.2	—	—	—	(BG)/(IR)	26.5
PLAPON60XO	1.42	0.15	0.17	Oil	Yes	(U)/BG/IR ***	26.5

Screw: W20.32 x 0.706 (0.8" x 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.


\*\*\*IR excitation might cause out of focus slightly.

**Universal Plan Semi Apochromat/Plan Semi Apochromat**


**UPLFLN, PLFLN series**

The UPLFLN and PLFLN objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. With their high S/N ratio, high resolution and contrast imaging, they are especially effective in brightfield and Nomarski DIC observations. For quality and performance, they offer an unbeatable solution to every kind of digital imaging need.


**UPLFLN4X**




**UPLFLN10X2**




**UPLFLN20X**



**UPLFLN40X**



**UPLFLN40XO**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
UPLFLN4X	0.13	17	—	—	—	(U)/(BG)/(IR)	26.5
UPLFLN10X2	0.30	10	—	—	—	U/BG/IR ***	26.5
UPLFLN20X	0.50	2.1	0.17	—	Yes	U/BG/IR ***	26.5
UPLFLN40X	0.75	0.51	0.17	—	Yes	U/BG/IR ***	26.5
UPLFLN40XO	1.30	0.2	0.17	Oil	Yes	U/BG/IR ***	26.5

Screw: W20.32 x 0.706 (0.8" x 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.


\*\*\*IR excitation might cause out of focus slightly.

## Universal Plan Semi Apochromat/Plan Semi Apochromat


### UPLFLN, PLFLN series

The UPLFLN and PLFLN objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. With their high S/N ratio, high resolution and contrast imaging, they are especially effective in brightfield and Nomarski DIC observations. For quality and performance, they offer an unbeatable solution to every kind of digital imaging need.


**UPLFLN60X**




**UPLFLN60XOI**




**UPLFLN100XO2**



**UPLFLN100XOI2**



**PLFLN100X**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
UPLFLN60X	0.90	0.2	0.11-0.23	—	Yes	U/BG/IR ***	26.5
UPLFLN60XOI	1.25-0.65	0.12	0.17	Oil	Yes	U/BG/IR ***	26.5
UPLFLN100XO2	1.30	0.2	0.17	Oil	Yes	U/BG/IR ***	26.5
UPLFLN100XOI2	1.3-0.6	0.2	0.17	Oil	Yes	U/BG/IR ***	26.5
PLFLN100X	0.95	0.2	0.14-0.2	—	Yes	BG	26.5

Screw: W20.32 × 0.706 (0.8" × 1/36")


\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.


\*\*\*IR excitation might cause out of focus slightly.

The PLN standard objectives provide excellent field flatness with fluorescence, darkfield and brightfield observation in transmitted light. These objectives are suited to clinical laboratory and examination work.


**PLN2X**




**PLN4X**



**PLN10X**



**PLN20X**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
PLN2X	0.06	5.8	—	—	—	(BG)	22.0
PLN4X	0.10	18.5	—	—	—	(BG)	22.0
PLN10X	0.25	10.6	—	—	—	(BG)	22.0
PLN20X	0.40	1.2	0.17	—	Yes	(BG)	22.0


Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.

The PLN standard objectives provide excellent field flatness with fluorescence, darkfield and brightfield observation in transmitted light. These objectives are suited to clinical laboratory and examination work.


**PLN40X**



**PLN50XOI**



**PLN100XO**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
PLN40X	0.65	0.6	0.17	—	Yes	(BG)	22.0
PLN50XOI	0.9-0.5	0.2	—	Oil	Yes	(BG)	22.0
PLN100XO	1.25	0.15	—	Oil	Yes	(BG)	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.

**Universal Plan Semi Apochromat for Phase Contrast**


**UPLFLN-PH series**

The UPLFLN-PH objectives are especially effective in phase-contrast observations with their high S/N ratio, high resolution and contrast imaging. These objectives display flat images from high transmission factors up to the near-infrared region of the spectrum.


**UPLFLN4XPH**




**UPLFLN10X2PH**




**UPLFLN20XPH**




**UPLFLN40XPH**



**UPLFLN60XOIPH**



**UPLFLN100XO2PH**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
UPLFLN4XPH	0.13	17	—	—	—	U/BG/IR ***	26.5
UPLFLN10X2PH	0.30	10	—	—	—	U/BG/IR ***	26.5
UPLFLN20XPH	0.50	2.1	0.17	—	Yes	U/BG/IR ***	26.5
UPLFLN40XPH	0.75	0.51	0.17	—	Yes	U/BG/IR ***	26.5
UPLFLN60XOIPH	1.25-0.65	0.12	0.17	Oil	Yes	U/BG/IR ***	26.5
UPLFLN100XO2PH	1.30	0.2	0.17	Oil	Yes	U/BG/IR ***	26.5

Screw: W20.32 × 0.706 (0.8" × 1/36")


\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.


\*\*\*IR excitation might cause out of focus slightly.

The PLN-PH standard objectives provide excellent field flatness with phase-contrast observation in transmitted light.


**PLN10XPH**




**PLN20XPH**



**PLN40XPH**



**PLN100XOPH**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
PLN10XPH	0.25	10.6	—	—	—	(BG)	22.0
PLN20XPH	0.40	1.2	0.17	—	Yes	(BG)	22.0
PLN40XPH	0.65	0.6	0.17	—	Yes	(BG)	22.0
PLN100XOPH	1.25	0.15	—	Oil	Yes	(BG)	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.

## Universal Plan Semi Apochromat for Polarizing


### UPLFLN-P series

The UPLFLN-P universal objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. These objectives reduce internal strain to a minimum and are proper for polarizing, Nomarski DIC, brightfield and fluorescence microscopy.


**UPLFLN4XP**




**UPLFLN10XP**




**UPLFLN20XP**



**UPLFLN40XP**



**UPLFLN100XOP**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
UPLFLN4XP	0.13	17	—	—	—	(U)/(BG)	26.5
UPLFLN10XP	0.30	10	—	—	—	U/BG	26.5
UPLFLN20XP	0.50	2.1	0.17	—	Yes	U/BG	26.5
UPLFLN40XP	0.75	0.51	0.17	—	Yes	U/BG	26.5
UPLFLN100XOP	1.30	0.2	0.17	Oil	Yes	U/BG	26.5

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.

**Achromat for Polarizing**


**PLN-P, ACHN-P series**

The PLN-P and ACHN-P cost-efficient objectives enable transmitted polarized light observations.


**PLN4XP**




**ACHN10XP**




**ACHN20XP**



**ACHN40XP**



**ACHN100XOP**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
PLN4XP	0.10	18.5	—	—	—	(BG)	22.0
ACHN10XP	0.25	6	—	—	—	(BG)	22.0
ACHN20XP	0.40	3	0.17	—	—	(BG)	22.0
ACHN40XP	0.65	0.45	0.17	—	Yes	(BG)	22.0
ACHN100XOP	1.25	0.13	—	Oil	Yes	(BG)	22.0

Screw: W20.32 x 0.706 (0.8" x 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.

**Plan Achromat (ND)**


**PLN-CY, PLFLN-CY series**

The PLN-CY and PLFLN-CY standard objectives provide excellent field flatness. These objectives, equipped with ND filters, provide the same level of brightness even if the magnification is changed, thereby removing the need for brightness adjustment.


**PLN2XCY**




**PLN4XCY**




**PLN10XCY**



**PLN20XCY**



**PLFLN10XCY**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
PLN2XCY	0.06	5.8	—	—	—	—	22.0
PLN4XCY	0.10	18.5	—	—	—	—	22.0
PLN10XCY	0.25	10.6	—	—	—	(BG)	22.0
PLN20XCY	0.40	1.2	0.17	—	—	(BG)	22.0
PLFLN10XCY	0.30	9.9	—	—	—	BG	26.5

Screw: W20.32 x 0.706 (0.8" x 1/36")


\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.

## Long Working Distance Universal Plan Semi Apochromat

### LUCPLFLN series

The LUCPLFLN long working distance, universal objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. These objectives are dedicated to tissue culture observations through bottles and dishes, offering high contrast and resolution in brightfield, DIC and fluorescence observations.



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
LUCPLFLN20X	0.45	7.8-6.6	0-2	—	—	U/BG/IR ***	22.0
LUCPLFLN40X	0.60	4-2.7	0-2	—	—	U/BG/IR ***	22.0
LUCPLFLN60X	0.70	2.2-1.5	0.1-1.3	—	—	U/BG/IR ***	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR : Near IR. "( )" indicates that fluorescence might be slightly darker.


\*\*\*IR excitation might cause out of focus slightly.

**Long Working Distance Universal Plan Semi Apochromat for Relief Contrast**


**CPLFLN-RC, LUCPLFLN-RC series**

The CPLFLN-RC and LUCPLFLN-RC long working distance, universal objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. These objectives are designed for the observation of living cells, including oocytes. Plastic vessels can be used with these objectives for relief-contrast observations.


**CPLFLN10XRC**



**LUCPLFLN20XRC**



**LUCPLFLN40XRC**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
CPLFLN10XRC	0.30	9 ***	—	—	—	BG	22.0
LUCPLFLN20XRC	0.45	7.8-6.6	0-2	—	—	U/BG/IR ****	22.0
LUCPLFLN40XRC	0.60	4.2-3.0	0-2	—	—	U/BG/IR ****	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.

\*\*\*Defined by 1 mm bottom thickness of Plastic Container, 0.5 mm bottom thickness of Glass Heat Plate (depends on the shape of container).


\*\*\*\*IR excitation might cause out of focus slightly.

**Long Working Distance Universal Plan Semi Apochromat for Phase Contrast**


**CPLFLN-PH, LUCPLFLN-PH series**

The CPLFLN-PH and LUCPLFLN-PH long working distance, universal objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. These objectives are exclusively designed for culture specimens and are optimized to produce exquisite phase-contrast images, regardless of the thickness and material of the vessel.


**CPLFLN10XPH**




**LUCPLFLN20XPH**



**LUCPLFLN40XPH**



**LUCPLFLN60XPH**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
CPLFLN10XPH	0.30	9.5 ***	1	—	—	BG	22.0
LUCPLFLN20XPH	0.45	7.8-6.6	0-2	—	—	U/BG/IR ****	22.0
LUCPLFLN40XPH	0.60	4.2-3.0	0-2	—	—	U/BG/IR ****	22.0
LUCPLFLN60XPH	0.70	2.2-1.5	0.1-1.3	—	—	U/BG/IR ****	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.

\*\*\*Defined by 1 mm bottom thickness of Plastic Container.


\*\*\*\*IR excitation might cause out of focus slightly.

The CPLN-PH and LCACHN-PH standard objectives provide excellent field flatness with phase-contrast observation in transmitted light.


**CPLN10XPH**



**LCACHN20XPH**



**LCACHN40XPH**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
CPLN10XPH	0.25	10 ***	—	—	—	BG	22.0
LCACHN20XPH	0.40	3.2 ***	1	—	—	BG	22.0
LCACHN40XPH	0.55	2.2 ***	1	—	—	BG	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")


\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.


\*\*\*Defined by 1 mm bottom thickness of Plastic Container.

The CPLN-RC and LCACHN-RC standard objectives provide excellent field flatness with relief-contrast observation in transmitted light. These objectives are designed for the observation of living cells, including oocytes, in plastic vessels.


**CPLN10XRC**



**LCACHN20XRC**



**LCACHN40XRC**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
CPLN10XRC	0.25	9.7 ***	—	—	—	BG	22.0
LCACHN20XRC	0.40	2.8 ***	1.5	—	—	BG	22.0
LCACHN40XRC	0.55	1.9 ***	1.5	—	—	BG	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.


\*\*\*Defined by 1 mm bottom thickness of Plastic Container, 0.5 mm bottom thickness of Glass Heat Plate (depends on the shape of container).

## No Cover Water Immersion for Fixed Stage Upright Microscope


### UMPLFLN-W, LUMPLFLN-W series

The UMPLFLN, LUMPLFLN-W and LUMFLN series of long working distance, water immersion objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. These objectives achieve exquisite DIC and fluorescence from the visible range to infrared and are ideal for fluorescence imaging of brain tissue, as well as other tissue and specimens.


**UMPLFLN10XW**




**UMPLFLN20XW**



**LUMPLFLN40XW**



**LUMPLFLN60XW**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
UMPLFLN10XW	0.30	3.50	0	Water	—	U/BG	26.5
UMPLFLN20XW	0.50	3.50	0	Water	—	U/BG	26.5
LUMPLFLN40XW	0.8	3	0	Water	—	U/BG	26.5
LUMPLFLN60XW	1.00	2	0	Water	—	U/BG	26.5

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.

**No Cover Water Immersion for Fixed Stage Upright Microscope**

**XLUMPLFLN20XW**

The XLUMPLFLN-W high NA, long working distance objectives display flat images from high transmission factors up to the near-infrared region of the spectrum. These objectives achieve exquisite DIC and fluorescence from the visible range to infrared. These objectives allow the measurement of cell membrane electric potential as the design of the objectives provides easy access to patch clamp electrodes.

**XLUMPLFLN20XW \*\*\***



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
XLUMPLFLN20XW ***	1.00	2	0	Water	—	U/BG/IR ****	22.0

Screw: W25 × 0.75

\*Defined with cover glass thickness shown as “CG” in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. “( )” indicates that fluorescence might be slightly darker.

\*\*\*Special dedicated nosepiece needed (WI-SNPXLU).


\*\*\*\*IR excitation might cause out of focus slightly.

**Universal Apochromat**


**UAPON 340 series**

The UAPON340 objectives feature high transmission of 340nm wavelength light, providing maximum performance in fluorescence microscopes through UV excitation.


**UAPON20XW340**



**UAPON40XO340-2**



**UAPON40XW340**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
UAPON20XW340	0.70	0.35	0.17	Water	Yes	U/BG	22.0
UAPON40XO340-2	1.35	0.1	0.17	Oil	Yes	U/BG	22.0
UAPON40XW340	1.15	0.25	0.13-0.25	Water	Yes	U/BG	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.


\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "( )" indicates that fluorescence might be slightly darker.

**TIRF Objectives**


**APON, UAPON series**

These objectives feature the highest NA to create an evanescent wave field for high-contrast TIRF images.


**APON60XOTIRF**




**APON100XHOTIRF \*\*\***



**UAPON100XOTIRF**



**UAPON150XOTIRF**



unit: mm

Objectives (Magnification)	NA	W.D. (mm) *	Cover Glass Thickness (mm)	Immersion	Spring	Fluorescence **	FN
APON60XOTIRF	1.49	0.1	0.13-0.19	Oil	—	(U)/BG	22.0
APON100XHOTIRF ***	1.70	0.08	0.15	Special Oil	—	BG	22.0
UAPON100XOTIRF	1.49	0.1	0.13-0.19	Oil	—	U/BG	22.0
UAPON150XOTIRF	1.45	0.08	0.13-0.19	Oil	—	U/BG	22.0

Screw: W20.32 × 0.706 (0.8" × 1/36")

\*Defined with cover glass thickness shown as "CG" in the above drawings.

\*\*U: UV excitation, BG: Visible Excitation, IR: Near IR. "()" indicates that fluorescence might be slightly darker.

\*\*\*Special Dedicated Cover Glass and Immersion Oil needed.




# Microscope Frame

## BX53M: Upright Transmitted & Reflected Light Microscope Frame

### BX53MTRF-S

Designed with modularity in mind, the BX3M series provide versatility for a wide variety of material science and industrial applications. The frames are outfitted with ESD capability to protect electronic samples.

BX53MTRF-S



Weight: 7.6 kg Unit: mm




# Microscope Frame

## BX53M: Upright Reflected Light Microscope Frame

### BX53MRF-S

Designed with modularity in mind, the BX3M series provide versatility for a wide variety of material science and industrial applications. The frames are outfitted with ESD capability to protect electronic samples.

BX53MRF-S



Weight: 7.4 kg Unit: mm




# Microscope Frame

## BX3: Automated Transmitted Light Microscope Frame

### BX63F

Fully motorized system allows automation of complex multidimensional experiments, The precise motorized Z-drive, High stability due to fixed stage design.

BX63F



Weight: 14.1 kg Unit: mm


# Microscope Frame

## BX3: Semi-Motorized Fluorescence Transmitted Light Microscope Frame

### BX53F

The entire optical path of the BX53 is designed for optimal fluorescence imaging, using UIS2 components that set new standards in precision and clarity. Modular concept allows motorization of individual components.

BX53F



Weight: 8.3 kg Unit: mm


# Microscope Frame

## BX3: Manual System Transmitted Light Microscope Frame

### BX43F

The BX43 microscope offers a wide range of features, exquisite optical performance and is the ideal platform for digital imaging. This flexible microscope offers various contrast methods and prominent optics combined with true-color LED illumination for remarkable color rendering.

BX43F



Weight: 9.1 kg Unit: mm




# Microscope Frame

## BX3: Transmitted Ergonomic Microscope Frame

### BX46F

BX46 features unrivalled ergonomic design with low-position fixed stage and nosepiece focus that increase working comfort.

BX46F






Weight: 9.8 kg Unit: mm

# Microscope Frame

## BX2: Upright Motorized Transmitted/Reflected Frame

BX61TRF

BX61TRF




Weight: 11.4 kg Unit: mm

# Microscope System BXFM

## BXFM Frame

### BXFM-F

Widely used system that allows use in combination with fiber illumination, motorized revolving nosepiece and telan lens unit. Can easily be integrated into other equipment. Attach to the equipment by rear bolt mounting screw or pillar mounting hole.




Weight: 1.9kg Unit: mm

# Microscope System BXFM

## BXFM System Configuration Example1


### BXFM-F + BXFM-ILH + BXFM-ILHSPU

Accommodates the reflected light brightfield/darkfield and fluorescence illuminators.



Weight: 3.2kg Unit: mm

### BXFM Combination Sample BXFM-F + BXFM-ILH + BXFM-ILHSPU + TR30-2 + BX-RLA2 + U-LH100L-3



Weight: 8.2kg (exclude objective) Unit: mm


\*For installation dimensions, refer to those for the BXFM-F.

# Microscope System BXFM

## BXFM System Configuration Example2


### BXFM-F + BXFM-ILHS

Compact focusing unit suitable for building into existing equipment.



Weight: 2.4kg Unit: mm

### BXFM-S Combination Sample BXFM-F + BXFM-ILHS + TR30-2 + U-KMAS + U-LH100L-3



\*For installation dimensions, refer to those for the BXFM-F.

Weight: 5.5 kg (exclude objective) Unit: mm


# Microscope System BXFM

## Stands for BXFM

A wide variety of stands are available to suit different applications and purposes.

### SZ2-STU2

Universal Stand Type 2




\*The rotation angle of the horizontal arm can restrict to 90 degrees with stopper.

#### Major Specifications

Item	Specifications
1 Diameter of Focusing Arm or Fixing Section of Tube	ø32 mm
2 Vertical Pole Diameter	ø40 mm
3 Horizontal Poles Diameters	ø25 mm (both upper and lower poles)
4 Stroke	Horizontal: 234 mm, Vertical: 205 mm
5 Movement Range	Horizontal: 541 (435 + 106) mm max. (vertical pole — BXFM-S optical axis)
6 Maximum Specimen Weight	Forward: 10 kg (within 90-degree area) Transverse Direction: 6 kg Backward Direction: 7 kg (at maximum stroke)
7 Weight	30 kg

### U-ST


Compact Stand



Weight: 1.8 kg

### SZ-STL

Large Stand



Weight: 5 kg


Unit: mm

# Illumination Units

## Reflected Light Illuminator for BX53M

### BX3M-RLA-S


Reflected Light Illuminator for BF/DF for BX53M



Weight: 3.0 kg

### BX3M-KMA-S

BF Reflected LED Light Illuminator for BX53M



Weight: 2.4 kg


Unit: mm

# Illumination Units

## Coded Reflected Light Illuminator for BX53M Frame

### BX3M-RLAS-S


Coded Reflected Light Illuminator for BF/DF for BX53M



Weight: 3.6 kg

### BX3M-URAS-S

Coded Universal Reflected Illuminator for BX53M



Weight: 3.2 kg

Unit: mm


# Illumination Units

## Reflected Illuminator for BX3 Series

### BX3-URA

Universal Reflected Illuminator for BX3 series

A total of eight fluorescence mirror units can be attached for comfortable multi-color fluorescence observations.




Weight: 3.8 kg

### BX3-RFAS

Coded Fluorescence Illuminator for BX3 series

Eight fluorescence mirror units with coded function.



Weight: 3.9 kg

Unit: mm

# Illumination Units

## Reflected Light Illuminator for BX2 series


### BX-RLA2

Reflected Light Illuminator for BF/DF

ND filters are linked when exchanging between brightfield and darkfield.

#### Accessories

Unit Name	Description	Weight (g)
U-25LBD	LBD Filter Slider	20
U-25IF550	IF550 Filter Slider	20
U-25ND6	ND Filter	20
U-25ND25	ND Filter	20
U-25FR	Frost Filter Slider	20
U-25L42	UV-cut Filter	20
U-PO3	Polarizer Slider for Reflected Light	71
U-POTP3	Polarizer Slider for Reflected Light with Tint Plate	71
U-AN360-3	360° Rotatable Analyzer Slider	79
U-AN	Analyzer Slider for Reflected Light	50
U-DICR	DIC Slider for Reflected Light	130
U-DICRH	DIC Slider for Reflected Light (Resolution type)	130
U-DICRHC	DIC Slider for Reflected Light (Contrast type)	130



Weight: 3.4 kg

### BX-URA2


Universal Reflected Light Illuminator

Suitable for observations ranging from brightfield to fluorescence.

Six mirror units can be attached to this reflected light illuminator simultaneously.

#### Accessories

Unit Name	Description	Weight (g)
U-25LBD	LBD Filter Slider	20
U-25IF550	IF550 Filter Slider	20
U-25ND6	ND Filter	20
U-25ND25	ND Filter	20
U-25FR	Frost Filter Slider	20
U-25L42	UV-cut Filter	20
U-PO3	Polarizer Slider for Reflected Light	71
U-POTP3	Polarizer Slider for Reflected Light with Tint Plate	71
U-AN360-3	360° Rotatable Analyzer Slider	79
U-AN	Analyzer Slider for Reflected Light	50
U-DICR	DIC Slider for Reflected Light	130
U-DICRH	DIC Slider for Reflected Light (Resolution type)	130
U-DICRHC	DIC Slider for Reflected Light (Contrast type)	130
U-MBF3	Mirror Unit for Reflected Brightfield	80
U-MDF3*	Mirror Unit for Reflected Darkfield	80
U-MDIC3	Mirror Unit for Reflected DIC	80
U-MBFL3	Mirror Unit for Reflected Brightfield, for High Intensity Light Source	80
U-MWUS3	Fluorescence Mirror Unit for Reflected (U excitation)	80
U-MWBS3	Fluorescence Mirror Unit for Reflected (B excitation)	80
U-MWGS3	Fluorescence Mirror Unit for Reflected (G excitation)	80




Weight: 3.8 kg

Unit: mm

# Illumination Units


Mounting Dimensions of Illuminator (BX3M-RLA-S, BX3M-RLAS-S, BX3M-URAS-S, BX3M-KMA-S, BX3-RFAS, BX3-URA)



Fix illuminator using four M5 screws and projection for fastening.

Unit: mm

Mounting Dimensions of Illuminators (BX-RLA2 and BX-URA2)



Fix illuminator using four M5 screws and projection for fastening.

Unit: mm

# Illumination Units

## Compact Reflected Light Illuminator for BF


### U-KMAS

Very compact reflected light illuminator with reduced depth.

### U-KMAS

#### Accessories

Unit Name	Description	Weight (g)
U-25LBD	LBD Filter Slider	20
U-25IF550	IF550 Filter Slider	20
U-25ND6	ND Filter	20
U-25ND25	ND Filter	20
U-25FR	Frost Filter Slider	20
U-25L42	UV-cut Filter	20
U-PO3	Polarizer Slider for Reflected Light	71
U-POTP3	Polarizer Slider for Reflected Light with Tint Plate	71
U-AN360-3	360° Rotatable Analyzer Slider	79
U-AN	Analyzer Slider for Reflected Light	50
U-DICR	DIC Slider for Reflected Light	130
U-DICRH	DIC Slider for Reflected Light (Resolution type)	130
U-DICRHC	DIC Slider for Reflected Light (Contrast type)	130



Weight: 1.2 kg


# Light Source Units

## LED Lamp Housing for BX53M

The BX3M utilizes a high-intensity white LED light source for both reflected and transmitted light. High-intensity light support various observation modes such as BF, DF, DIC and Polarizing.

### BX3M-LEDR


LED Lamp Housing for BX3M Reflected



Weight: 0.5 kg

### BX3M-LEDT


LED Lamp Housing for BX3M Transmitted



Weight: 0.5 kg

### BX3M-PSLED

Power Supply for LED Lamp House (refuired with BXFM)



Weight: 0.36 kg

Unit: mm


# Light Source Units

## LED Lamp Housing for BX3 Series

The BX3M utilizes a high-intensity white LED light source for both reflected and transmitted light. High-intensity light support various observation modes such as BF, DF, DIC and Polarizing.

### U-LHLEDC

LED Lamp Housing for BX3M Transmitted




Weight: 0.5 kg

### U-LEDPS

External Power Supply for U-LHLEDC

The BX3 utilizes a high color reproducible LED illumination that support a constant color temperature for consistent and support color reproduction of brightfield stains.



Weight: 0.78 kg

Unit: mm

# Light Source Units


## Lamp Housings

Various different lamp housings are available, for use with different light sources: choose to suit the intended purpose.

\*If you use the units in your production line, please consult your nearest Olympus representative in your region about the use conditions beforehand.

### U-LH75XEAP0

75W Xenon Apo Lamp Housing




\*Power supply unit (U-RX-T) and power cable (UYCP) are necessary for 75 W xenon lamp housing. These items are sold separately.

Cable Length 2,000 mm Accepted Lamp: USH-1030L Weight: 2.7 kg

### U-LH100HGAP0

100W Mercury Apo Lamp Housing



### U-LH100HG

100W Mercury Lamp Housing

\*Power supply unit (U-RFL-T) and power cable (UYCP) are necessary for 100 W mercury lamp housings. These items are sold separately.


Cable Length 2,000 mm Accepted Lamp: USH-1030L Weight: 2.7 kg

### U-RFL-T

Power Supply Unit for Mercury Lamp

### U-RX-T

Power Supply Unit for Xenon Lamp



Weight: Approximately 3 kg

Unit: mm

# Light Source Units


## Halogen Illumination

For the 100 W halogen lamp, the external power supply TH4-100/200 with an intensity adjustment switch and an ON/OFF switch, both are located close to the operator's hand, are provided.

\*If you use the units in your production line, please consult your nearest Olympus representative in your region about the use conditions beforehand.

*Illumination devices for microscope have suggested lifetimes. Periodic inspections are required. Please visit our web site for details.*

### U-LH100-3/U-LH100IR/U-LH100L-3 100 W Halogen Lamp Housings




Cable Length U-LH100-3: 290 mm  
 U-LH100IR: 290 mm  
 U-LH100L-3: 800 mm

Accepted Lamp: 12V100WHAL (high intensity lamp)  
 12V100WHAL-L (long life lamp)

\*External power supply (TH4-100 or TH4-200) and power cable (UYCP) are necessary for 100 W halogen lamp housings. These items are sold separately.


Weight: 880 g

### TH4-100/200 External Power Supply




Weight: 2.2 kg

### TH4-HS Hand Switch



Weight: 140 g

### U-RMT Extension Cord



Weight: 200 g

Unit: mm

# Light Source Units

## Halogen Fiber Illumination Accessories


All Olympus reflected light illuminators can be used with fiber illumination.

\*If you use the units in your production line, please consult your nearest Olympus representative in your region about the use conditions beforehand.

*Illumination devices for microscope have suggested lifetimes. Periodic inspections are required. Please visit our web site for details.*

### U-RCV


DF Converter for BX-URA2



Weight: 315 g

### U-LGAD


Fiber Adapter for Reflected Light Observation



Light Guide Mount Hole  $\phi 12$  Weight: 390 g

### LG-PS2\*

Light Source




\*The types of model varies by country in use.

Weight: 1.6 kg

### LG-SF

Light Guide



Weight: 210 g

Unit: mm


# Light Source Units

## Lamp Housing Accessory

Two lamp housings can be attached simultaneously.

\*If you use the units in your production line, please consult your nearest Olympus representative in your region about the use conditions beforehand.

**U-DULHA**  
Double Lamp House Adapter




Weight: 1.2 kg



# Condenser Units


## Condenser

**U-AC2-7**  
Abbe Condenser





Weight: 174 g

**U-LWCD**  
Long Working Distance Condenser





Weight: 380 g

**WI-OB CD**  
Long Working Distance Oblique Condenser



Weight: 530 g

**CH3-CD**  
Brightfield Condenser



Weight: 60 g

Unit: mm


# Observation Tubes

## Super Widefield Trinocular Observation Tubes

Trinocular observation tubes with super widefield of view. Compatible with FN 26.5.


### U-SWTR-3

Super Widefield Trinocular Tube



### MX-SWETTR

Super Widefield Erect Image Tilting Trinocular Tube



Unit: mm

Name	FN	Inclination Angle (degree)	Interpupillary Distance (mm)	Light Path Selector (eyepiece/video port)	Observation Image	Weight (g)
U-SWTR-3	26.5	24	50-76	100/0, 20/80, 0/100	Inverted	2300
MX-SWETTR	26.5	0-42	50-76	100/0, 0/100	Erect	4200

\*Length marked with an asterisk (\*) may vary according to interpupillary distance. The distance for figure shown is 62 mm.

# Observation Tubes

## Widefield Trinocular Observation Tubes


Trinocular observation tubes with widefield of view. Compatible with FN 22.

### U-TR30-2

Widefield Trinocular Tube


### U-TR30IR

Widefield Trinocular Tube for IR



### U-ETR-4

Widefield Erect Image Trinocular Tube



Unit: mm

Name	FN	Inclination Angle (degree)	Interpupillary Distance (mm)	Light Path Selector (eyepiece/video port)	Observation Image	Weight (g)
U-TR30-2	22	30	50-76	100/0, 20/80, 0/100	Inverted	1600
U-TR30IR	22	30	50-76	100/0, 0/100	Inverted	1600
U-ETR-4	22	30	50-76	100/0, 0/100	Erect	1900

\*Length marked with an asterisk (\*) may vary according to interpupillary distance. The distance for figure shown is 62 mm.

## Single Port Tube with Lens


When the visual observation is not needed and only video observation is required, a single port tube with a built-in telan lens can be attached directly to the video port.

### U-TLU

Single Port Tube with Lens

### U-TLUIR

Single Port Tube with Lens for IR



• For attachable video camera adapters, refer to video camera adapters system diagram page (pages 2-12 - 2-13).


Weight: 350 g Unit: mm

# Observation Tubes

## Tilting Binocular Tube


### U-TTBI

Telescopic Tilting Binocular Tube




### U-TTR-2

Tilting Trinocular Tube



### U-TBI-3

Tilting Binocular Tube



Unit: mm


Name	FN	Inclination Angle (degree)	Interpupillary Distance (mm)	Light Path Selector (eyepiece/video port)	Observation Image	Weight (g)
U-TTBI	22	0-25	50-76	NA	Inverted	3800
U-TTR-2	22	5-35	50-76	100/0, 0/100, 50/50	Inverted	3200
U-TBI-3	22	5-35	50-76	NA	Inverted	1300

\*In the case of combination with CX-RFA-2, FN reduces to 18.


# Observation Tubes

## Binocular Tube

**U-CBI30-2**  
Binocular Tube



**U-CTR30-2**  
Trinocular Tube



Unit: mm

Name	FN	Inclination Angle (degree)	Interpupillary Distance (mm)	Light Path Selector (eyepiece/video port)	Observation Image	Weight (g)
U-CBI30-2	20, 18(*)	30	48-75	NA	Inverted	800
U-CTR30-2	20, 18(*)	30	48-75	NA	Inverted	900

\*In the case of combination with CX-RFA-2, FN reduces to 18.

# Intermediate Tubes & Accessories


## Intermediate Tubes

Various accessories for various observation need.

### U-CA

Magnification Changer

Provides 1X, 1.2X, 1.6X and 2X intermediate magnifications.



Weight: 1.3 kg

### U-ECA


Magnification Changer 2C

Provides 1X and 2X intermediate magnifications.

### U-ECA1.6X

Magnification Changer 1.6X

Provides 1X and 1.6X intermediate magnifications.




Weight: 1.3 kg

### U-TRU

Trinocular Intermediate Attachment

Intermediate attachment which divides the light path, allowing attachment of both digital and video cameras.



BI:PT=100:0/20:80

Weight: 1.3 kg

Unit: mm

# Intermediate Tubes & Accessories


## Intermediate Tubes

Various accessories for various observation need.

### U-DP

Dual Port

Use this intermediate tube to divide the light path.



Weight: 1 kg

Light Path Selector by Mirror Unit



Transmitted Side Port: Side Port = 100:0




Transmitted Side Port: Side Port = 70:30 (with use of U-MBF3)

### U-DP1XC

Dual Port 1X

Combine with U-DP to obtain a 1X image.




Weight: 500 g

### U-EPA2

Eyepoint Adjuster


Raises eyepoint by 30 mm.



Weight: Approximately 500 g

### U-APT

Arrow Pointer



Weight: 1.2 kg

Unit: mm


# Intermediate Tubes & Accessories

## Dual port tube with C mounts

### U-DPCAD

Dual camera port allows the attachment of two cameras.

U-DPCAD




Weight: 0.9 kg Unit: mm

# Eyepieces


## Eyepieces

Eyepieces for UIS2 optical system.


**WHN10X**  
Widefield Eyepiece




**WHN10X-H**  
**CROSSWHN10X**  
Widefield Eyepieces



**WH15X**  
Widefield Eyepiece



**SWH10X-H**  
**MICROSSWH10X**  
**CROSSSWH10X**  
Super Widefield Eyepieces



unit: mm

Name	FN	Diopter Adjustment Range (1/m)	Micrometer Diameter (mm)	Weight (g)	Remarks
WHN10X	22	—	24	90	
WHN10X-H	22	-8 — +5	24	170	with Adjustable Diopter
CROSSWHN10X	22	-8 — +5	—	170	with Cross Lines and Adjustable Diopter
WH15X	14	—	24	90	
SWH10X-H	26.5	-8 — +2	—	210	with Adjustable Diopter
MICROSSWH10X	26.5	-8 — +2	—	210	with Micrometer and Adjustable Diopter
CROSSSWH10X	26.5	-8 — +2	—	210	with Cross Lines and Adjustable Diopter

\*EP=eyepoint


# Revolving Nosepieces

## Revolving Nosepieces for BF Objectives

Choose from following 6 types. For motorized nosepieces, refer to motorized unit page.

### U-5RE-2


Quintuple Revolving Nosepiece for BF/DF with Slider Slot for DIC



Weight: 520 g

### U-D6RE


Septuple Revolving Nosepiece with Slider Slot for DIC



Weight: 800 g

### U-D6RE-ESD


Sextuple Revolving Nosepiece with Slider Slot for DIC with ESD Treatment



Weight: 0.7 kg

### U-D7RE


Sextuple Revolving Nosepiece with Slider Slot for DIC



Weight: 980 g

### U-P4RE


Centerable Quadruple Revolving Nosepiece with Slider Slot for DIC



Weight: 1 kg

### U-P6RE

Centerable Sextuple Revolving Nosepiece with Slider Slot for DIC



Weight: 1 kg

Insert the DIC dummy when not using the DIC slider. unit: mm


# Revolving Nosepieces

## Revolving Nosepieces for BF/DF Objectives

Choose from following 3 types. Use of adapter to mount BF objectives (BD-M-AD) enables attachment of brightfield objectives. For motorized nosepieces, refer to motorized unit page.

### U-5BDRE


Quintuple Revolving Nosepiece for BF/DF



Weight: 800 g

### U-D5BDRE


Quintuple Revolving Nosepiece for BF/DF with Slider Slot for DIC



Weight: 800 g

### U-P5BDRE

Centerable Quintuple Revolving Nosepiece with Slider Slot for DIC




Weight: 1 kg

Insert the DIC dummy when not using the DIC slider.

### U-D6BDRE

Sextuple Revolving Nosepiece for BF/DF with Slider Slot for DIC

### BD-M-AD Adapter to Mount BF Objectives



Weight: 10 kg


unit: mm

# Revolving Nosepieces

## Coded Sextuple Revolving Nosepiece

### U-D6RES


Coded Sextuple Revolving Nosepiece  
with Slider Slot for DIC



Weight: 0.7 kg

### U-D6BDRES-S


Coded Sextuple BF/DF Revolving  
Nosepiece with Slider Slot for DIC



Weight: 0.9 kg

### U-D5BDRES-ESD


Coded Quintuple BF/DF Revolving  
Nosepiece with Slider Slot for DIC with  
ESD Treatment



Weight: 0.75 kg

### U-D7RES

Coded Septuple Nosepiece with Slider Slot for DIC



unit: mm


# Video Camera Adapters

## C-mount Video Camera Ports

Allows direct attachment of a C-mount video camera. Five types are provided: 1X, 0.63X, 0.5X, 0.35X and 0.25X. All models feature a focus adjustment function.

### U-TV0.25XC


C-mount Video Port with 0.25X Lens



Weight: 1.2 kg

### U-TV0.35XC-2


C-mount Video Port with 0.35X Lens



Weight: 100 g

### U-TV0.5XC-3


C-mount Video Port with 0.5X Lens



Weight: 200 g

### U-TV0.63XC


C-mount Video Port with 0.63X Lens




Weight: 430 g

### U-TV1XC

C-mount Video Port with 1X Lens



Weight: 300 g



Video Camera Adapter (projection lens)	Projection Magnifications	Projection Area (FN)		
		2/3" CCD	1/2" CCD	1/3" CCD
U-TV1X-2	1X	11	8	6
U-TV1XC	1X	11	8	6
U-TV0.63XC	0.63X	17.5	12.7	9.5
U-TV0.5XC-3	0.5X	22	16	12
U-TV0.35XC-2	0.35X	—	22	17.1
U-TV0.25XC	0.25X	—	—	24

$$\text{Practical Field of View (mm)} = \frac{\text{Projection Area (FN)}}{\text{Objective Magnifications}}$$

Focus the video camera adapter to prevent defocusing the eyepiece image and defocusing by magnification switching. Generally, the video camera adapter is focused by switching to a low magnification after focusing at a high magnification objective.


Unit: mm

# Video Camera Adapters

## Video Camera Mount Adapters


Allows attachment to video cameras with C, Bayonet, Sony and F mounts. Use with the U-TV1X-2. Focus by amount of screwing into U-TV1X-2.

**U-CMAD3**  
C-mount Adapter




Weight: 165 g

**U-BMAD**  
Bayonet Mount Adapter




Weight: 80 g

**U-SMAD**  
Sony Mount Adapter




Weight: 90 g

**U-TMAD**  
T Mount Adapter



Weight: 165 g

**U-FMT**  
F/T Mount Adapter \*




\*It must be combined with U-TMAD.

Weight: 80 g

## Video Camera Port

This port can be attached directly to the trinocular observation tube as well as to the single port tube with lens.

**U-TV1X-2**  
Video Port 1X



Weight: 80 g

Unit: mm


# Motorized Units

## Motorized Universal Reflected Illuminator for BX2 Series

Various motorized units, appropriate for automation of equipment, are available.

### BX-RLAA + U-D6REMC + U-LH100-3

Motorized BF/DF Reflected Light Illuminator + Motorized Nomarski DIC Sextuple Revolving Nosepiece + 100 W Halogen Lamp Housing  
 Enables motorized exchange of objectives, selection between brightfield and darkfield observations as well as aperture diaphragm closing/opening. The BX-UCB control unit has an RS232C connector, allowing control via a PC. For method of attaching illuminator, refer to Mounting Dimensions of BX-RLA2.




Illuminator Cable Length: 1.800 mm Weight: 5.5 kg (exclude objective)

### BX-RFAA

#### Motorized Universal Reflected Light Illuminator

Reflected light fluorescence illuminator with simultaneous attachment of six mirror units. Incorporates motorized mirror unit changeover and shutter.



Illuminator Cable Length: 1.800 mm Weight: 4.3 kg

Unit: mm

# Motorized Units

## Motorized Units


Various motorized units, appropriate for automation of equipment, are available.

### U-P5BDREMC

Motorized Centerable Quintuple BD Revolving Nosepiece with Slider Slot for DIC

### U-D6BDREMC

Motorized Sextuple BD Revolving Nosepiece with Slider Slot for DIC



Weight: 80 g

### U-D5BDREMC


Motorized Quintuple BD Revolving Nosepiece with Slider Slot for DIC

### U-D6REMC

Motorized Sextuple Revolving Nosepiece with Slider Slot for DIC

### U-P5REMC

Motorized Centerable Quintuple Revolving Nosepiece with Slider Slot for DIC



Unit: mm

# Motorized Units


## Control Box for BX2 Series

Various motorized units, appropriate for automation of equipment, are available.

### BX-UCB

Control Unit

Motorized units including motorized illuminator and auto focus unit can be totally controlled from BX-UCB.




\*Extension cord U-RMT (1700 mm) should be used to connect the lamp housing (U-LH100-3) to the BX-UCB.

Weight: 1.0 kg

### U-HSTR2

Hand Switch




Cable Length 2000 mm Weight: 370 g

### BX-REMCB

Control Box for Motorized Nosepiece and BF/DF Illuminator

BX-RLAA and U-D5BDREMC/U-D6REMC/U-P5REMC can be controlled from U-HSTR2, or direct from the computer keyboard via an RS232C connector.


\*BX-RFAA and U-D5BDREM/U-D6REM combination not applicable.



Weight: kg

### U-ACAD4515

AC Adapter for BX-REMCB



Weight: kg

Unit: mm

# Motorized Units

## Motorized Units for BX2 Series


Various motorized units, appropriate for automation of equipment, are available.

### U-AFA2M-VIS

Active Auto Focus Unit

Featuring an AF laser light source in wavelength 785 nm. The multiple-spots sensor enables the high-speed and stable focusing of specimens with variable height differences

\*Class 1 laser product



\*Consult your Olympus dealer about the motorized focus.

Weight: 2.6 kg

Cable


AFA2M-CBL2M Weight: 0.23 kg Length: 2 m

AFA2M-CBL3M Weight: 0.36 kg Length: 3 m

### U-FWR

Motorized Reflected Filter Wheel

Accomplish maximum 6 filter position exchange.



Weight: 1.0 kg

Unit: mm

# Motorized Units


## Motorized Modular Microscope

A motorized microscope unit for integration with your equipment. Motorized operations such as revolving nosepiece up/down, objective switching, aperture diaphragm open/close, and brightfield/darkfield switching are accomplished with this component. Several microscopic operations are totally controlled from an external unit by combining this component with an auto focus unit.

### BXFMA-F

#### Motorized Illumination with Power Focus

This is the configuration combined with Active Auto Focus Unit U-AFA2M-VIS, Single Port Tube with Lens U-TLU, a lamp housing, a motorized nosepiece and objectives.




\*Consult your Olympus dealer about the mounting dimensions.

Weight: 13 kg (BXFMA-F frame 7.6 kg)

### U-FH


#### Focus Adjustment Knob Unit



Weight: 760 g

### U-IFFH

#### Focus Adjustment Knob Interface



Weight: 1450 g


Unit: mm

# Motorized Units

## Motorized Units for BX3 Series

The flexibility of the motorized fluorescence illuminator accommodates multi-color 'stained specimens. The 8-position mirror units permit quick changeover of fluorescence colors.

**BX3-RFAA**  
Motorized Fluorescence Illuminator




Weight: 4.2 kg

**BX3-UCD8A**  
Motorized Universal Condenser


The motorized universal condenser integrates a variety of optical elements to accommodate transmitted light techniques including brightfield, DIC and phase contrast observation.

\*Optical Elements are required. Please consult your nearest Olympus representative in your region about the Optical Elements Combination.



Weight: 1.6 kg

**U-D7REA**  
Motorized Septuple Revolving Nosepiece with Slider Slot for DIC



Weight: 1.5 kg


Unit: mm

# Motorized Units

## Control Box for BX3 Series

### BX3-CBH

Control Box for Fully-motorized Function  
Control Box connected to PC via IEEE1394




Weight: 4.2 kg

### U-CBM

Control Box M for Motorized Function  
Control box connected to PC via RS232

\*Please consult your nearest Olympus representative in your region about the detailed system chart.




Weight: 0.82 kg

### U-CBS

Control Box for Coded Function

Control Box connected to PC via RS232 (attached onto Microscope Frame)

\*Please consult your nearest Olympus representative in your region about the detailed system chart.



Weight: 0.5 kg

Unit: mm

# Motorized Units




## Control Box for BX53M/BXFM

### BX3M-CB

Control Box for BXFM

Control Box connected to PC via RS232/USB2.0 (attached onto Microscope Frame)

\*Please consult your nearest Olympus representative in your region about the detailed system configuration and combination.






Weight: 0.49 kg

### BX3M-CBFM

Control Box for BX53M Frame

Control Box connected to PC via RS232/USB2.0

\*Please consult your nearest Olympus representative in your region about the detailed system configuration and combination.



Weight: 0.67 kg

Unit: mm

# Optical Terminology

## 1. FN and Practical Field of View

The field number (FN) is referred to as the diaphragm size of eyepiece in mm unit which defines the image area of specimen. The diaphragm diameter actually seen through eyepiece is known as the practical field of view (F.O.V.) which is determined by the formula:


$$\text{F.O.V.} = \frac{\text{Eyepiece FN}}{\text{Objective Magnification}} \text{ (mm)}$$

## 2. W.D.

The distance between the front edge of the objective and the specimen surface (with the surface of the cover glass in case of the cover glass objective) when the specimen is focused.

## 3. Parfocal Distance

It is the distance between the objective mounting plane and the specimen. In UIS2/UIS objectives, the parfocal distance is designed at 45 mm.



For parfocal distance of the LCPLFLN-LCD series objectives, refer to the appropriate objective page.

## 4. Relationship Between the Objective's Focal Length and Magnifications

Indicated magnifications of UIS2/UIS objectives are the values when the focal length of the tube lens is 180 mm.

$$M_{(ob)} = \frac{\text{Focal Length of Tube Lens}}{f}$$

$M_{(ob)}$ : Objective Magnification  
 f: Objective's Focal Length

## 5. Total Magnification

### 5.1 Observation Through Eyepiece (binocular observation)

$$M_{(bino)} = M_{(ob)} \times M_{(oc)}$$

$M_{(bino)}$ : Total magnification for binocular observation  
 $M_{(ob)}$ : Objective magnification  
 $M_{(oc)}$ : Eyepiece magnification

### 5.2 Video Monitor Observation

#### ● Total Magnification for Video Monitor

$$M_{(video\ monitor)} = M_{(ob)} \times M_{(video\ camera\ adapter)} \times \text{Monitor Magnification}^*$$

$M_{(video\ monitor)}$ : Total Magnification on the Video Monitor  
 $M_{(ob)}$ : Objective Magnification  
 $M_{(video\ camera\ adapter)}$ : Projected Magnification for Video Camera Adapter Including Photo Eyepiece  
 (refer to Figure 1)

\*Refer to Figure 3 for "Monitor Magnification"

#### ● Practical Field of View for Video Monitor Observation

$$\text{Practical Field of View for Video Monitor Observation} = \frac{\text{Image Device Size}^*}{M_{(ob)} \times M_{(video\ camera\ adapter)}}$$

$M_{(ob)}$ : Objective Magnification  
 $M_{(video\ camera\ adapter)}$ : Projected Magnification for Video Camera Adapter Including Photo Eyepiece  
 (refer to Figure 1 for projected magnifications)

\*Refer to Figure 2 for Image Device Size

Figure 1 Video Camera Adapter and Projection Magnifications

Video Camera Adapter (projection lens)	Projection Magnifications
U-TV1X-1 + Video Camera Mount Adapters	1X
U-TV1XC	1X
U-TV0.63XC	0.63X
U-TV0.5XC-3	0.5X
U-TV0.35XC-2	0.35X
U-TV0.25XC	0.25X

Figure 2 Imaging Device Size

Camera Format	Diagonal	Horizontal	Vertical
1/3"	6.0 mm	4.8 mm	3.6 mm
1/2"	8.0 mm	6.4 mm	4.8 mm
2/3"	11.0 mm	8.8 mm	6.6 mm

The above table is for standard image device sizes. Check your device size for precise calculation.

Figure 3 Imaging Device Size and Monitor Magnifications

Camera Format	Monitor Size (diagonal)				
	10"	15"	17"	19"	21"
1/3"	42.3X	63.5X	72.0X	80.4X	88.9X
1/2"	31.8X	47.6X	54.0X	60.3X	66.7X
2/3"	23.1X	34.6X	39.3X	43.9X	48.5X

#### Example

What is total magnifications for video monitor when objective is 50X, video camera adapter U-TV0.5XC, 2/3" video camera and 21" monitor are used ?

# Optical Terminology

•Total magnification on the video monitor:

$m_{(ob)}=50X$ ,  $M_{(video\ camera\ adapter)}$  is 0.5X from Figure 1 and monitor magnification is 48.5x from Figure 3.

$M_{(monitor\ observation)}=M_{(ob)} \times M_{(video\ camera\ adapter)} \times \text{monitor magnification}$   
 $=50 \times 0.5 \times 48.5=1213 \times$

•Practical field of view for video observation (horizontal side):

$M_{(ob)}=50X$ ,  $M_{(video\ camera\ adapter)}$  is 0.5X from Figure 1 and horizontal side of 2/3" imaging device is 8.8mm from Figure 2

$$\begin{aligned} \text{Practical Field of View} &= \frac{\text{Image Device Size}}{M_{(ob)} \times M_{(video\ camera\ adapter)}} \\ \text{for Video Observation} &= \frac{8.8\text{ (mm)}}{50 \times 0.5} = 52\ \mu\text{m} \end{aligned}$$

## 6. NA

The numerical aperture is a key factor to the performance of objective (resolving power, focal depth and brightness). The NA is determined by the following formula:


$$NA = n \times \sin\theta$$

$n$ =Refraction Rate of the Medium between Specimen and Objectives. (Air:  $n=1$ , oil:  $n=1.515$ )

$\theta$ : Angle Which is Made by the Optical Axis and Refraction of the Light Farthest from the Center of the Lens.

The visual field brightness (B) of the microscope is determined by the following formula in relation to the objective magnification (M). The larger the NA and the lower the objective magnification, brightness will increase in the factor of the second power.

$$B \propto \frac{NA^2}{M^2}$$



## 7. Resolving Power

The resolving power of an objective is measured by its ability to differentiate two lines or points in an object. The greater the resolving power, the smaller the minimum distance between two lines or points that can still be distinguished. The larger the NA, the higher the resolving power.

## ●Resolving Power Formula

The following formula is generally used for determining resolution.

$$\epsilon = 0.61 \times \frac{\lambda}{NA} \quad (\text{Reyleigh formula})$$

$\lambda$ : Wavelength or Radiation in Use

( $\lambda=0.55\ \mu\text{m}$  is used for visible light.)

NA: Objective NA

Example

MPLFLN100X (NA=0.90),  $\lambda=0.55\ \mu\text{m}$

$$\epsilon = 0.61 \times \frac{\lambda}{NA} = \frac{0.3355}{NA} = \frac{0.3355}{0.90} = 0.37\ \mu\text{m}$$

## 8. Focal Depth of Microscope

The focal depth refers to the depth of the specimen layer which is in sharp focus at the same time, even if the distance between the objective and the specimen plane is changed when observing and shooting the specimen plane by microscope. As human eyes are individually different in the ability of their focus adjustment, each person's perception of the focal depth varies.

At present, the Berek formula is generally used, because it gives a focal depth value that often coincides with that obtained through experiments.

### Focal Depth formula

#### ●Visual Observation (Berek formula)

$$\pm \text{D.O.F.} = \frac{\omega \times 250,000}{NA \times M} + \frac{\lambda}{2(NA)^2} \quad (\mu\text{m})$$

D.O.F.: Depth Of Focus

$\omega$ : Resolving Power of Eyes 0.0014

(when optical angle is 0.5 degrees)

M: Total Magnification

(objective magnification x eyepiece magnification)

$$\rightarrow \pm \text{D.O.F.} = \frac{350}{NA \times M} + \frac{0.275}{NA^2} \quad (\lambda=0.55\ \mu\text{m})$$

This indicates that the focal depth becomes smaller as the numerical aperture becomes larger.

Example

With MPLFLN100X (NA=0.90), WHN10X:

$$\pm \text{D.O.F.} = \frac{350}{0.90 \times 1,000} + \frac{0.275}{0.81} = 0.39 + 0.34 = 0.73\ \mu\text{m}$$

#### ●Video Camera

In the case of a video camera, the focal depth will vary according to number of pixels of CCD, optical magnification, and numerical aperture. The above-mentioned formula is used as a rough guide only.


## 9. Aberrations

A difference between an ideal image and an actual image that passes through an optical system is called an “aberration.”

### 9.1 Requirements for Ideal Image Formation

The following three requirements must be satisfied to form an image with no aberration, or an ideal image.

- (i) All the light rays coming from a single point and passing through an image formation optical system converge on a single point.
- (ii) Image points, which correspond to object points on the same plane perpendicular to the optical axis, are present on the same plane.
- (iii) The planar shape of an object and the planar shape of an image that are on the same plane perpendicular to the optical axis have a similarity relation.




In an actual optical system, however, it is very difficult to strictly meet the requirements for ideal image formation and this causes “aberrations” that interfere with image forming performance.

### 9.2 Classification of Aberrations

Aberrations that interfere with image forming performance are classified as shown below in Figure 9-2.

Seidel's Aberration = “Expansion of a Point Image” + “Curvature of Image Plane” + “Deformation”




Types (1) to (3) correspond to “expansion of a point image” that goes against requirement (i) for ideal image formation in Figure 9-1. Type (4) corresponds to “curvature of image plane” that goes against requirement (ii) in Figure 9-1.

Type (5) corresponds to “deformation” that goes against requirement (iii) in Figure 9-1.

Types (6) and (7) correspond to “color blur” of images caused by characteristics of glass materials used for the optical system. “Expansion of a point image” can also be expressed by “wavefront aberration” that regards the light as “waves” and takes account of the phase to include the influence of diffraction.

#### (1) Spherical Aberration


When light rays coming out of an axial object point enter a lens, the light rays with a larger numerical aperture (NA) are subjected to stronger refraction power and cross the optical axis in positions with larger differences from the ideal image formation position. The aberration caused this way by different image forming positions due to differences in NA of axial light rays is called “spherical aberration.” (“Spherical aberration” is proportional to the cube of NA)



It is said that objectives with larger NA have better resolution but worsen spherical aberration. Our advanced design and manufacturing techniques have realized good optical performance even with large numerical aperture.

#### (2) Coma Aberration


Even though spherical aberration is compensated to be very small, there are cases where light rays coming out of an off-axis object point are not condensed to a single point on the image plane but generate asymmetric blur just like a comet leaving traces. This is called coma aberration.



# Optical Terminology

## (3) Astigmatism

Even though a lens is compensated for spherical aberration and coma aberration, there are cases where an image of an off-axis object point is not focused to a single point but separated to a concentric line image and a radial line image. This is called "astigmatism." When astigmatism is present, a point image blurs vertically and horizontally, before and after the focus position.




## (4) Field Curvature

An image plane of an object on a plane perpendicular to an optical axis does not always become a plane perpendicular to the optical axis, but it generally becomes a curved plane. This symptom is called "field curvature."

When field curvature is present, the image is more displaced as it becomes closer to the periphery of the visual field. Therefore, when the center of an image is brought into focus, blur occurs in the peripheral areas of the image. To bring the entire image, including the periphery, into clear focus, it is necessary to adequately compensate for this type of aberration.

## (5) Distortion

When there is no similar relation between a planar shape on an object and a shape on the image plane, this is called "distortion." When distortion is present, a square image appears in a shape of a barrel or pin-cushion as shown in Figure 9-6.



The microscope optical system contains some distortion. When distortion is present, it can bring erroneous results of shape measurements. When a microscope is used for precision measurements, pay close attention to this aberration, for example, by providing it with an aberration compensation function.

## (6) Chromatic Aberration


Glasses used for optical systems have different refractive indexes depending on the wavelength. This causes differences in focal length between wavelengths and generates displacement of image forming position. This phenomenon is called "chromatic aberration," which is sometimes subdivided into axial displacement on the optical axis, called "axial chromatic aberration" (or lateral chromatic aberration) and displacement on the image plane, called "chromatic aberration of magnitude."

Many special glass materials are used, e.g., for apochromats, to eliminate chromatic aberration in a wide range from violet light (g-rays with wavelength of 435 nm) to red light (c-rays with wavelength of 656 nm).


### 9.3 Wavefront Aberration

Since a long time ago, aberrations have been used in “geometric optics,” which considers light as “light rays.” Microscope optical systems are often used for observation of very small specimens at a wavelength level, and sometimes adopt “wave optics,” which regards light as “waves” and handles the phase information, taking account of the influence of diffraction.

In such a case, “wavefront aberration” is used for evaluation. As shown below, when requirements for ideal imaging are satisfied in a microscope optical system, the spherical wavefront (spherical waves) coming from a single point on an object (specimen) is converted to plane waves through an ideal objective. The plane waves are converted to spherical waves through an ideal tube lens, and condensed to a single point. The wavefront of these waves is called the “ideal wavefront.”




Based on the figure indicated for (1) spherical aberration, the behavior of the wavefront in an optical system that has an aberration is described below.




A difference (a degree of disagreement) between the ideal wavefront and the actual wavefront shown above is called “wavefront aberration.”

### 9.4 Strehl Ratio

When a point light source is observed with an aberration-free optical system and an aberrated optical system, the former concentrates the focal point to a point at the image formation position. In contrast, the latter fails to produce a focal point, instead causing a spread in the intensity distribution of the point image (this is known as “point spread”). The specific appearance of such a point image (i.e. point spread) is shown in Figure 9-9.



With the proportion of light concentrated in the image plane (intensity of light concentrated in the Airy disk) by an aberration-free optical system serving as 100%, the proportion of light concentrated by an aberrated optical system is known as the Strehl ratio. When graphed, the Strehl ratio reveals peaks in intensity as shown in Figure 9-10. The higher the S.R., the closer an optical system is to being aberration-free.



A Strehl ratio of 80% is typically called the diffraction limit, and lenses with a lower ratio lack the performance required to serve as an objective. A ratio of over 95% means that the lens' performance in general observations is comparable to that of an aplanatic lens (which is corrected for spherical aberrations and coma).

*Note: A laser interferometer is used for actual assessment of optical performance, so assessment is done at a single wavelength. Unless otherwise noted, Strehl ratio measurements are at the e-line (544 nm).*

- **OLYMPUS CORPORATION is ISO14001 certified.**
- **OLYMPUS CORPORATION is ISO9001 certified.**
- All company and product names are registered trademarks and/or trademarks of their respective owners.
- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.
- Illumination devices for microscope have suggested lifetimes. Periodic inspections are required.  
Please visit our web site for details.

<http://www.olympus-ims.com/en/component/>

**OLYMPUS**<sup>®</sup>

For enquiries - contact  
[www.olympus-ims.com/contact-us](http://www.olympus-ims.com/contact-us)

**OLYMPUS CORPORATION**

Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

**OLYMPUS EUROPA SE & CO. KG**

Wendenstrasse 14-18, 20097 Hamburg, Germany

**OLYMPUS SCIENTIFIC SOLUTIONS AMERICAS CORP.**

48 Woerd Avenue, Waltham, MA 02453, U.S.A.

**OLYMPUS MEDICAL SYSTEMS INDIA PRIVATE LIMITED.**

Ground Floor, Tower-C, SAS Tower, The Medicity Complex, Sector-38, Gurgaon 122001, Haryana, India

**OLYMPUS KOREA CO., LTD.**

8F Olympus Tower, 446 Bongeunsa-ro, Gangnam-gu, Seoul, 06153 Korea

**OLYMPUS AUSTRALIA PTY. LTD.**

3 Acacia Place, Notting Hill VIC 3168, Australia

**OLYMPUS SINGAPORE PTE LTD.**

491B River Valley Road, #12-01/04 Valley Point Office Tower, Singapore 248373

**OLYMPUS (CHINA) CO., LTD.**

10F, K. Wah Centre, 1010 Huai Hai Road (M), Xuhai District, Shanghai, 200031 P.R.C.